

Stretching the limits of impedance testing



1260

Impedance/gain-phase Analyzer

The 1260 Impedance/gain-phase Analyzer is - without doubt - the most powerful, accurate and flexible Frequency Response Analyzer available today.

In daily use by leading researchers wherever measurement integrity and experimental reliability are of paramount importance, 1260's solid reputation is frequently endorsed in published research papers in fields such as:-

- Corrosion studies
- · Battery research and fuel cells
- Solar cells
- LCDs
- Bio-materials
- Ceramics / composites
- Electronic component development
- · Civil engineering

Part of Solartron Analytical's extensive range of precision products designed to provide cost effective solutions for dc and ac analysis in electrochemical and materials research, 1260 offers an outstanding measurement specification for impedance spectroscopy:

Huge frequency range

Spanning $10\mu Hz$ to 32MHz with 0.015ppm resolution, 1260 provides excellent coverage for virtually all chemical and molecular mechanisms - all in a single instrument.

Unbeatable accuracy

With an accuracy of 0.1%, 0.1°, measurements can be made with complete confidence, and even the most subtle changes in sample behavior detected and quantized.

Noise free analysis

1260 uses Solartron Analytical's patented single-sine correlation technique, which inherently removes the noise and harmonic distortion which plagues lesser instruments.

- Frequency resolution: 1 in 65 million (0.015ppm)
- 0.1%, 0.1° accuracy unsurpassed by any similar instrument
- Resolution to 0.001dB, 0.01° capturing every detail
- Measures impedances $>100M\Omega$
- 2-, 3- and 4-terminal measurement configurations
- Polarization voltage up to ±40.95V
- Renowned ZPlot software package simplifies experiments and optimises throughput

Systems

When combined with other products from Solartron Analytical's range, including well-proven application software, 1260 can form the heart of an advanced electrochemical and materials measurement system, to provide superb accuracy, flexibility and reliability - even for the most complex research problems.

Impedance measurement

Virtually every liquid and solid is able to pass current when a voltage is applied to it. If a variable (ac) voltage is applied to the material, the ratio of voltage to current is known as the impedance. The measured impedance varies with the frequency of the applied voltage in a way that is related to the properties of the liquid or solid. This may be due to the physical structure of the material, to chemical processes within it or a combination of both.

The advantages of impedance measurement over other techniques include:-

- · Rapid aquisition of data
- Accurate, repeatable measurements
- Non-destructive
- Highly adaptable to a wide variety of different applications.
- Ability to differentiate effects due to electrodes, diffusion, mass/charge transfer by analysis over different frequency ranges
- Equivalent circuit/modelling techniques for detailed analysis of results



Please visit us at: www.valuetronics.com

1260 Impedance/gain-phase Analyzer Specification

Frequency

Generator Voltage mode Current mode 0 to 60mA rms ac Amplitude ≤10MHz 0 to 3V rms 0 to 1V rms 0 to 20mA rms >10MHz

Maximum ac resolution 5mV 100µA ±40.95V $\pm 100 mA$ dc bias range 10mV 200µA Maximum dc resolution Output impedance 50Ω±1% $>200k\Omega$ at <1kHz

error: ±100ppm, stability, 24hrs ±1°C: ±10ppm frequency (log or lin), ac/dc voltage, ac/dc current Sweep types

range: 10μHz to 32MHz, max resolution: 10μHz

Maximum voltage hi to lo: ±46V peak, lo to ground: ±0.4V peak Maximum current ±100mA peak

Impedance lo to ground: $100k\Omega$, <10nFConnection single BNC, floating shield Output disable contact closure or TTL logic 0

Input System Voltage (2x) Current

3 independent analyzers operating in parallel Ranges 30mV, 300mV, 3V 6μA, 60μA, 600μA, 6mA, 60mA Maximum resolution 1µV 200pA Full scale peak ±5V ±100mA

Inputs protected to ±46V ±250mA Connections single/differential BNC single BNC

Shields floating/grounded

Coupling dc or ac (-3dB at 1Hz) dc or ac (-3dB at Hz)

Input impedance Hi to shield 1Mohm, <35pF \geq 600µA range, 1Ω

10kohm, 330pF $<600\mu A$ range, 50Ω Shield to ground

Limits of error Ambient temperature 20±10°C, integration time >200ms.

Data valid for one year after calibration.

Results Variable frequency, ac amplitude, dc bias

Measured parameters voltage gain, phase, real, imaginary, Z, R, X, Y, G, B, V, I


group delay, C, L, Q, D

Power supply 90 to 126V, 198 to 252V, 48 to 65Hz Power consumption 230VA

432mm x 176mm x 573mm (17in x 6.93in x 22.56in) Dimensions (w x h x d)


Weight 18kg (40lbs)

Operating temp. range 0 to 50°C (32 to 122°F) Limit of error



Impedance Measurements


Applies for stimulation level of 1V for impedances $>50\Omega$ or 20mA for impedances <50 Ω



Solartron Analytical is a world leader in instrumentation and software for the characterization of materials and electrochemical systems using precision electrical measurement techniques.

These techniques find particular use in the fields of corrosion, battery and fuel cell research, dielectric analysis and electrochemistry. The product portfolio includes industry standard frequency response analyzers, potentiostats, electrochemical software (Zplot and CorrWare) and battery test equipment.

Arun Technology, an operating unit of Solartron Analytical, provides a range of metal analyzers using optical emission techniques for determining elemental content. The units in static laboratory or mobile format are used in foundries, steelworks, or scrapyards for metals analysis or material identification.



Windows is a trademark of Microsoft Corporation ZPlot and CorrWare are trademarks of Scribner Associates Inc.



Gain-phase

measurements

Applies to all ranges

at >10% full scale