

Agilent N77xx Series

N7744A 4-Channel Optical Multiport Power Meter (LXI) N7745A 8-Channel Optical Multiport Power Meter (LXI)

User's Guide



Notices

© Agilent Technologies, Inc. 2010

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Manual Part Number

N7744-90B01

Edition

Fourth edition, July 2010

Agilent Technologies Deutschland GmbH Herrenberger Str. 130 71034 Böblingen, Germany

Warranty

This Agilent Technologies instrument product is warranted against defects in material and workmanship for a period of one year from date of shipment. During the warranty period, Agilent will, at its option, either repair or replace products that prove to be defective.

For warranty service or repair, this product must be returned to a service facility designated by Agilent. Buyer shall prepay shipping charges to Agilent and Agilent shall pay shipping charges to return the product to Buyer. However, Buyer shall pay all shipping charges, duties, and taxes for products returned to Agilent from another country.

Agilent warrants that its software and firmware designated by Agilent for use with an instrument will execute its programming instructions when properly installed on that instrument. Agilent does not warrant that the operation of the instrument, software, or firmware will be uninterrupted or error free.

Limitation of Warranty

The foregoing warranty shall not apply to defects resulting from improper or inadequate maintenance by Buyer, Buyer-supplied software or interfacing, unauthorized modification or misuse, operation outside of the environmental specifications for the product, or improper site preparation or maintenance.

No other warranty is expressed or implied. Agilent Technologies specifically disclaims the implied warranties of Merchantability and Fitness for a Particular Purpose.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with the terms of such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

Exclusive Remedies

The remedies provided herein are Buyer's sole and exclusive remedies. Agilent Technologies shall not be liable for any direct, indirect, special, incidental, or consequential damages whether based on contract, tort, or any other legal theory.

Assistance

Product maintenance agreements and other customer assistance agreements are available for Agilent Technologies products. For any assistance contact your nearest Agilent Technologies Sales and Service Office.

Certification

Agilent Technologies Inc. certifies that this product met its published specifications at the time of shipment from the factory. Agilent Technologies further certifies that its calibration measurements are traceable to the United States National Institute of Standards and Technology, NIST (formerly the United States National Bureau of Standards, NBS) to the extent allowed by the Institutes's calibration facility, and to the calibration facilities of other International Standards Organization members.

ISO 9001 Certification

Produced to ISO 9001 international quality system standard as part of our objective of continually increasing customer satisfaction through improved process control.

Safety Notices

CAUTION

A **CAUTION** notice calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in damage to the product or loss of important data. Do not proceed beyond a **CAUTION** notice until the indicated conditions are fully understood and met.

WARNING

A WARNING notice calls attention to an operating procedure, practice, or the like that, if not correctly performed or adhered to, could result in personal injury or death. Do not proceed beyond a WARNING notice until the indicated conditions are fully understood and met.

Contents

1 Getting Started

```
Overview
          7
Agilent N7744A and N7745A Optical Power Meter 7
Agilent N775xA Optical Attenuator and Power Meter and N776xA Optical Attenuator
  8
Agilent N7711/4A Single and 4-Port Tunable Laser System Source
                                                               10
Initial Inspection 10
Safety Considerations 11
  Safety Considerations - Overview 11
  Safety Symbols
                   11
  General 12
  Operating Environment 13
  Line Power Requirements
                            13
  Initial Safety Information for Tunable Laser Modules 13
  Laser Safety Labels
                       14
  Input/Output Signals
                         15
  Line Power Connectors 15
  Instrument Markings
                        16
AC Line Power Supply Requirements 17
  AC Line Power Supply Requirements - Overview
                                                 17
  Line Power Requirements 17
  Line Power Cable 17
  Changing the Fuse 19
Operating and Storage Environment 19
  Operating and Storage Environment - Overview
                                                19
  Temperature
               19
  Humidity 20
  Altitude
            20
  Pollution Protection
                      20
  Instrument Cooling
                      20
```

Agilent N77xx Series

```
Input and Output Connectors 21
  Input and Output Connectors - Overview
                                         21
  Front Panels - Multiport Power Meters
                                        21
  Front Panels - Optical Attenuators 22
  Front Panels - Tunable Laser Source
                                     24
  Rear Panel 24
Making Optical Connections 25
  Making Optical Connections
                              25
  Making Quad-adapter Connections on the Agilent N7744/45A Multiport Power
  Meter 26
Optical Output
               27
Electrical Connectors
                     28
LAN Interface 28
  LAN Interface
                 28
  Selecting a LAN Network 29
  Private LAN Considerations 29
  Security 29
  Performance
                29
  Reliability 30
  IP Address Availability 30
  Site LAN Considerations 30
  Connecting the LAN Cable
                            30
GPIB Interface 33
  GPIB Interface 33
  Cables and Adapters 33
  Connector 33
  GPIB Logic Levels
                     34
USB Interface
              34
Powering Up the Instrument 35
  Powering Up the Instrument 35
  Front Panel Indicators 35
Connecting to the Instrument 36
  Connecting to the Instrument - Overview
                                         36
  Connecting over USB
                        36
  Finding the IP Address of an instrument 37
  Connecting over LAN
                        39
Claims and Repackaging 42
  Claims and Repackaging
                           42
  Return Shipments to Agilent Technologies 42
```

N77xx Viewer Common User Interface Functions 43 Connecting to an Instrument 43 Viewing Measurements 45 Saving and Recalling Configurations 46 Controlling the Instrument Refresh 47

2 **Power Measurement**

Power Measurement 49 Controlling the Optical Multiport Power Meter 49 How to Measure Power 53 The Power Value 53 How to Set the Power Unit 54 How to Set the Calibration Offset 55 How to Set the Reference Level **55** How to Set the Wavelength 57 How to Remove Electrical Offsets (Zero) 57 How to Choose the Range Mode 58 How to Set the Range 59 How to Set the Averaging Time **60** How to Set Automatic Gain for Averaging Modulated Signals 60 How to Choose the MinMax Mode 60 How to Use External Triggers 62

3 Power Meter Applications

Power Meter Applications65The Logging Application65Continuous Logging Using Data Buffers66The Stability Application67

4 Specifications and Regulations Compliance

Definition of Terms 69 Literature 88 Technical Specifications - Optical Multiport Power Meter 89 Declarations of Conformity - Optical Multiport Power Meter 90 Regulations Information 91 EMC Canada 92 Acoustic Noise Information Germany 92

Agilent N77xx Series

Remove all doubt 92 Remove all doubt 92 Agilent E-mail Updates 92 Agilent Direct 93 Agilent Open 93 LXI 93 www.agilent.com 94

5 Accessories

Accessories - Overview 95 Instrument and Options - Agilent N7744A and N7745A 95 Cables 96

6 Maintenance and Troubleshooting

Maintenance and Troubleshooting - Overview 97 Cleaning Instructions 97 Cleaning Instructions 97 Cleaning Instructions for this Instrument 98 Safety Precautions 98 Why is it important to clean optical devices? 98 What do I need for proper cleaning? 99 Preserving Connectors 104 Cleaning Instrument Housings 104 Which Cleaning Procedure should I use? 105 How to clean connectors 105 How to clean connector adapters 107

How to clean connectors 105 How to clean connector adapters 107 How to clean bare fiber adapters 107 How to clean optical glass plates 108 Additional Cleaning Information 108 Other Cleaning Hints 109 Firmware Upgrades 109 Error Messages 110 SYST:ERR? 110 Data Sanitation 110

Memory Inventory 110 Clearing User Data 110

Index 113



Overview

This chapter helps you start using your N77-Series instrument.

- "Agilent N7744A and N7745A Optical Power Meter" on page 7
- "Agilent N775xA Optical Attenuator and Power Meter and N776xA Optical Attenuator" on page 8
- "Initial Inspection" on page 10
- "Safety Considerations Overview" on page 11
- "AC Line Power Supply Requirements Overview" on page 17
- "Operating and Storage Environment Overview" on page 19
- "Input and Output Connectors Overview" on page 21
- "Making Optical Connections" on page 25
- "Electrical Connectors" on page 28
- "LAN Interface" on page 28
- "GPIB Interface" on page 33
- "USB Interface" on page 34
- "Powering Up the Instrument" on page 35
- "Connecting to the Instrument Overview" on page 36
- "Claims and Repackaging" on page 42

Agilent N7744A and N7745A Optical Power Meter

The Agilent N7744A and Agilent N7745A are LXI compliant optical power meters with 4 or 8 power-sensor ports respectively. They address productivity challenges by increasing throughput and operational efficiency using industry-leading device connectivity, high-speed



Agilent Technologies

acquisition of measurement data and fast data transfer for postprocessing.

The Agilent N7744A and Agilent N7745A are designed to characterize optical multiport components, such as multiplexers, PON splitters, wavelength selective switches (WSS) and ROADMs, in all product phases from R&D to manufacturing.

The one rack unit (1U half-rack high) instrument can be operated alone or with other units, by PC connected over LAN, USB, or GPIB.

The unprecedented N7740A fiber connectivity concept consists of a quadruple adapter (quad-adapter) with a quick-locking mechanism. Up to four fiber connectors from the device under test (DUT) can be connected in a comfortable working position to one quad-adapter. The push-on locking mechanism makes it easy and secure to attach the quad-adapter to and remove the quad-adapter from the instrument in seconds. With multiple quad-adapters, the measurement task can be fully separated from the DUT connection task. This means a further increase in test throughput. Where this is not an issue, traditional single connector connectivity is supported as well.

These power meters are highly compatible with setups and programs developed using the Agilent 816x modular power meters. They are also easy to synchronize with Agilent tunable laser sources, to make spectral measurements. The functionality for time-dependent measurement has been extended with sampling at up to 1 MHz and data logging of up to 1 million samples for each port.

Agilent N775xA Optical Attenuator and Power Meter and N776xA Optical Attenuator

The Agilent N775xA and N776xA are compact multichannel attenuators and power meters. They are a new class of remote-controlled fiber optic instruments. All attenuators feature power control functionality to set the output power level of the attenuator.

The Agilent N7751A and Agilent N7752A are optical attenuators for single mode fiber with one and two single-mode fiber attenuator channels respectively and two extra optical power meter channels, for convenient power measurements and setup calibration. The Agilent N7761A, Agilent N7762A and Agilent N7764A are single-mode fiber optical attenuators for single mode fiber, with one, two and four attenuator channels respectively.

All N775xA and N776xA attenuators address productivity challenges by increasing throughput and operational efficiency.

The Agilent N7766A and the Agilent N7768A are optical attenuators for multimode fibers with two or four attenuator channels repectively. Both multimode attenuators are available with $50\mu m$ and $62.5\mu m$ multimode fiber.

The N7766A and N7768A are 2 and 4 channel multimode fiber attenuators that are optimized for quick power setting and preservation of the input modal power distribution.

The attenuators are designed to attenuate and control the power of light in optical fibers. They allow you to accurately and repeatably set the attenuation factor or power level manually using the Viewer software, or using the programming interface.

The attenuators include power control functionality that allow you to set the output power level of the attenuator. When power control mode is enabled, the module automatically corrects for power changes at the input so that the output power level you set is maintained.




Figure 1 Optical Attenuators with Power Control

The one rack unit (1U half-rack high) instrument can be operated alone or with other units, controlled over LAN, USB, or GPIB.

These attenuators and power meters are fully compatible with setups and programs developed using the Agilent 8157x modular attenuators and basic 8163x optical power meter commands, not including logging functionality.

Modal Fidelity for Multimode Fiber Systems

Signals in multimode fibers are distributed over a range of mode groups that can have different loss and delay in a link. For dependable multimode transceiver testing, the instrument used to set the power level should not change this modal distribution. The bulk-optic filter and collimated beam path of Agilent multimode attenuators are designed to assure homogeneous attenuation of all input modes.

Agilent N77xx Series

Agilent N7711/4A Single and 4-Port Tunable Laser System Source

The Agilent N7711A and N7714A tunable lasers are single-port and fourport sources, available with C-band or L-band wavelength coverage. The narrow linewidth and offset grid fine-tuning capability of the N7711A and N7714A make them ideal sources for realistic loading of the latest transmission systems.

To owners of Agilent's proven Lightwave Measurement System the 81950A tunable laser source module offers the same features as the N7711A. The 81950A plugs into the 8163B and 8164B mainframes. For additional information about the 81950A please refer to the Compact Tunable Laser data sheet, publication no. 5988-8518EN.

All models can reach any wavelength point within their specified wavelength range just like all other Agilent tunable lasers. In this mode, code compatibility with existing test setups based on Agilent's range of full-size and compact tunable lasers is a great asset. In system loading applications, it may be preferable to grid-tune the lasers like system transmitters, simply by changing the channel index. The channel grid is adjustable to standard ITU-T grid spacing like 50 GHz, and to arbitrary grids. Likewise, the zero frequency (base channel) of the chosen grid is adjustable. A 12 GHz fine-tuning range allows de-tuning the frequency.

Targeted for high test throughput, lowest cost-per-channel, and narrow footprint, all members of Agilent's 77-Series optical test instruments are built on a common platform and a common PC-based user interface. A complete set of control interfaces including LAN, USB2.0 and GPIB simplifies integration with manufacturing control systems. Code compatible to Agilent's Lightwave Measurement System modules, the new instrument generation can serve as plug-in replacements in existing test solutions. The N77 viewer software allows to control all features of the N7711/4A tunable laser sources.

Initial Inspection

Inspect the shipping container for damage. If there is damage to the container or cushioning, keep them until you have checked the contents of the shipment for completeness and verified the instrument both mechanically and electrically. The Function Tests give a procedure for checking the operation of the instrument. If the contents are incomplete, mechanical damage or defect is apparent, or if an instrument does not pass the operator's checks, notify the nearest Agilent Technologies Sales/Service Office.

NOTE

To avoid hazardous electrical shock, do not perform electrical tests when there are signs of shipping damage to any portion of the outer enclosure (covers, panels, and so on).

Safety Considerations


Safety Considerations - Overview

The following general safety precautions must be observed during all phases of operation, service, and repair of this instrument. Failure to comply with these precautions or with specific warnings elsewhere in this manual violates safety standards of design, manufacture, and intended use of the instrument. Agilent Technologies Inc. assumes no liability for the customer's failure to comply with these requirements. This product has been designed and tested in accordance with IEC Publication 61010-1, Safety Requirements for Electrical Equipment for Measurement, Control and Laboratory, and has been supplied in a safe condition. The instruction documentation contains information and warnings that must be followed by the user to ensure safe operation and to maintain the product in a safe condition.

Safety Symbols

	product. Do not proceed beyond a caution sign until the indicated conditions are fully understood and met.
W A R N I N G	The <i>warning</i> sign denotes a hazard. It calls attention to a procedure which, if not correctly performed or adhered to, could result in injury or loss of life. Do not proceed beyond a warning sign until the indicated conditions are fully understood

Agilent N77xx Series



General

This is a Safety Class 1 instrument (provided with a protective earth terminal) and has been manufactured and tested according to international safety standards.

Before operation, you should review the instrument and manual for safety markings and instructions. You must follow these to ensure safe operation and to maintain the instrument in safe condition.

Some circuits are powered whenever the instrument is connected to the AC power source. To disconnect from the line power, disconnect the power cord either at the rear power inlet or at the AC line power source (receptacle). One of these must always be accessible. If the instrument is in a cabinet, it must be disconnected from the line power by the system's line power switch.

WARNING

To avoid hazardous electrical shock, do not perform electrical tests when there are signs of shipping damage to any portion of the outer enclosure (covers, panels, and so on).

Operating Environment

WARNING

The instrument is not designed for outdoor use. To prevent potential fire or shock hazard, do not expose the instrument to rain or other excessive moisture.

Line Power Requirements

CAUTION

The instrument complies with installation category II and can operate from the single-phase AC power source that supplies between 100 V and 240 V at a frequency in the range 50 to 60 Hz. The maximum voltage fluctuation is 10% of the nominal supply voltage. The maximum power consumption is 60 VA with all options installed.

Initial Safety Information for Tunable Laser Modules

The laser sources specified by this user guide are classified according to IEC 60825-1 (2007) Laser Notice No. 50 dated 2007-June-24.

The laser sources comply with 21 CFR 1040.10 except for deviations pursuant to Laser Notice No. 50 dated 2001-July-26:

Table 2

	Agilent N7711A, N7714A, N7721A, N7722A	Agilent N7711A, N7714A, N7721A, N7722A	
	Options 110, 111, 210, 211, 222, 240	Options 101, 111, 201, 211, 222, 204	
Laser Type	ECL-Laser InGaAsP	ECL-Laser InGaAsP	
Wavelength range	1527nm-1566nm	1570nm-1609nm	
Max. CW output power*	50 mW	50 mW	
Beam waist diameter	9 µm	9 µm	

Agilent N77xx Series

	Agilent N7711A, N7714A, N7721A, N7722A	Agilent N7711A, N7714A, N7721A, N7722A	
	Options 110, 111, 210, 211, 222, 240	Options 101, 111, 201, 211, 222, 204	
Numerical aperture	0.1	0.1	
Laser Class according to IEC 60825-1 (2007)- Intl.	1 M	1 M	
Max. permissible CW output power**	163 mW	163 mW	

* Max. CW output power is defined as the highest possible optical power that the laser source can produce at its output connector.

** Max. permissible CW output power is the highest optical power that is permitted within the appropriate laser class.

Laser Safety Labels

Laser class 1M label




Figure 2 Class 1M Safety Label - N7711/14A, N7721/22A

A sheet of laser safety labels is included with the laser module as required. In order to meet the requirements of IEC 60825-1 we recommend that you stick the laser safety labels, in your language, onto a suitable location on the outside of the instrument where they are clearly visible to anyone using the instrument.

14

WARNING

Please pay attention to the following laser safety warning: Under no circumstances look into the end of an optical cable attached to the optical output when the device is operational. The laser radiation can seriously damage your eyesight. Do not enable the laser when there is no fiber attached to the optical output connector. The laser is enabled by pressing the 'active' button close to the optical output connector on the front panel of the module. The laser is on when the green LED on the front panel of the instrument is lit. The use of optical instruments with this product will increase eye hazard. The laser module has a built-in safety circuitry which will disable the optical output in the case of a fault condition Refer servicing only to qualified and authorized personnel.

Input/Output Signals

CAUTION

There is one input BNC connector: Trigger In. This is a TTL input. A maximum of 5 V can be applied as an external voltage to this input connector.

There is one output BNC connector: Trigger Out. This is a TTL output. Do not apply an external voltage to this connector.

Line Power Connectors

In accordance with international safety standards, the instrument has a three-wire power cable. When connected to an appropriate AC power receptacle, this cable earths the instrument cabinet. The type of power cable shipped with each instrument depends on the country of destination. Please refer to "Accessories - Overview" on page 95" for the part numbers of available power cables.

Agilent N77xx Series

WARNING


To avoid the possibility of injury or death, you must observe the following precautions before switching on the instrument.

- Insert the power cable plug only into a socket outlet provided with a protective earth contact. Do not negate this protective action by the using an extension cord without a protective conductor.
- Do not interrupt the protective earth connection intentionally.
- Do not remove protective covers. Operating personnel must not remove instrument covers. Component replacement and internal adjustments must be made only by qualified service personnel.
- Instruments that appear damaged or defective should be made inoperative and secured against unintended operation until they can be repaired by qualified service personnel.
- Defective, damaged, or malfunctioning laser sources must be returned to an Agilent Technologies Service Center.

Do not operate the instrument in the presence of flammable gases or fumes. Operation of any electrical instrument in such an environment constitutes a definite safety hazard.

Instrument Markings

	The instruction manual symbol. The product is marked with this warning symbol when it is necessary for the user to refer to the instructions in the manual.
Œ	The CE mark is a registered trademark of the European Community.
€₽°	The CSA mark is a registered trademark of the Canadian Standards Association.
C	The C-Tick mark is a registered trademark of the Australian Spectrum Management Agency.



AC Line Power Supply Requirements

AC Line Power Supply Requirements - Overview

This secton provides information on:

- Line Power Requirements
- Line Power Cable
- Changing the Fuse

Line Power Requirements

The instrument complies with installation category II and can operate from the single-phase AC power source that supplies between 100 V and 240 V at a frequency in the range 50 to 60 Hz. The maximum voltage fluctuation is 10% of the nominal supply voltage. The maximum power consumption is 60 VA with all options installed.

Line Power Cable

In accordance with international safety standards, the instrument has a three-wire power cable. When connected to an appropriate AC power receptacle, this cable earths the instrument cabinet. For the part number of the power cable for your country and instrument see the "Accessories - Overview" on page 95 chapter.


Agilent N77xx Series

CAUTION

Please note that the switch on the front panel of the instrument does not stop the flow of power to the instrument.

If you need to turn off the power, unplug the instrument at the mains or remove the power cable connector from the appliance coupler at the rear of the instrument. For this reason, the power cable connection should be easily accessible - allowing you to turn off the power quickly. If the instrument is in a cabinet, it must be disconnected from the line power by the system's line power switch.

The power switch allows you to switch between stand-by mode and power-on mode.





When the instrument is in stand-by mode, the Pwr LED is orange. When the instrument is powered-on, the Pwr LED is green.

WARNING

To avoid the possibility of injury or death, you must observe the following precautions before switching on the instrument.

- Insert the power cable plug only into a socket outlet provided with a protective earth contact. Do not negate this protective action by using an extension cord without a protective conductor.
- Do not interrupt the protective earth connection intentionally.

The AC power requirements are summarized on the rear panel of the instrument.




Figure 4 AC Power Requirement Markings

Changing the Fuse

CAUTION

There is no user replaceable fuse for the instrument. Changing the fuse should be carried out only by Agilent Technologies service personnel. If you need to get the fuse replaced refer to your nearest Agilent Technologies Sales/Service Office.

Operating and Storage Environment



The instrument should be protected from temperature extremes and changes in temperature that may cause condensation within it.

The operating temperature is from 5° C to $+40^{\circ}$ C.

The storage temperature is from -40 $^{\circ}\mathrm{C}$ to +70 $^{\circ}\mathrm{C}.$

Agilent N77xx Series

Humidity

The operating humidity is 15 to 95%, non-condensing.

Altitude

The maximum operating altitude is 2000 m.

Pollution Protection

The instruments are designed for pollution degree 2.

Instrument Cooling

The instrument has a cooling fan mounted internally.

Mount or position your instrument upright and horizontally, as shown in Figure 4 so that air can circulate through it freely.

Operating Position

When operating the instrument choose a location that provides at least 75 mm (3 inches) of clearance at the rear, and at least 25 mm (1 inch) of clearance at each side. Failure to provide adequate air clearance may result in excessive internal temperature, reducing instrument reliability. The instrument should not be operated when resting on its rear or side panels.




Figure 5 Correct Operating Position

(shown here for the Agilent N7744/5A, valid also for the N7751/2A and Agilent N7761/2/4A, N7766/68A and the N7711/4A)

Input and Output Connectors

Input and Output Connectors - Overview

This section provides information on:

- Front Panels Multiport Power Meters
- Front Panels Optical Attenuators
- Front Panels Tunable Laser Source
- Rear Panel

Front Panels - Multiport Power Meters

Agilent N77xx Series




Figure 6 Front panel of the Agilent N7744A 4-Channel Optical Multiport Power Meter




Figure 7 Front panel of the Agilent N7745A 8-Channel Optical Multiport Power Meter

Front Panels - Optical Attenuators




Figure 8 Front panel of the Agilent N7751A 1-Channel Optical Attenuator and 2-Channel Power Meter





Figure 9 Front panel of the Agilent N7752A 2-Channel Optical Attenuator and 2-Channel Power Meter








Figure 11 Front panel of the Agilent N7762A and N7766A 2-Channel Optical Attenuator




Figure 12 Front panel of the Agilent N7764A and N7768A 4-Channel Optical Attenuator

Agilent N77xx Series

Front Panels - Tunable Laser Source




Figure 13 Front panel of the Agilent N7711A Single-Port Tunable Laser System Source

Agiliant Technologies	NT714A Multiport Tunabis Lassr Source	
OBroor OBtace O OLAN LAN C Reset	Laser 1 On - O PMF PMF PMF PMF PMF PMF PMF On - O PMF PMF On - O On - O On - O On - O PMF On - O On	Laser 4 On -O PMP Cellend (D) Out
Lann		

Figure 14 Front panel of the Agilent N7714A Single-Port Tunable Laser System Source

Rear Panel




Figure 15 Rear panel of the Agilent N7744/5A Optical Multiport Power Meter, Agilent N7751/2A Optical Attenuator and Power Meter, and Agilent N7761/2/4/6/8A Optical Attenuator and Agilent N7711/4A Tunable Laser System Source

Making Optical Connections

Making Optical Connections

The optical attenuators are designed for use with FC connectors.

 $Attenuators \, ordered \, with \, option \, 021 \, have \, straight \, connector \, interfaces.$

Attenuators ordered with option 022 have angled connector interfaces.

The tunable lasers are designed for use with FC connectors.

Lasers ordered with option 071 have straight connector interfaces.

Lasers ordered with option 072 have agled connector interfaces.

The N7751/2A power meters also use FC connectors but are noncontact. They can be used with both angled and straight connectors.

The N7766A and N7768A multimode attenuators have straight connector interfaces. The fiber option determines the dimensions of the internal fiber core, such as 50 or $62.5\mu m$. Best performance is achieved by matching this to the fibers used in the setup.

CAUTION

If the connectors on your attenuator or lasers are angled, you can only use cables with angled connectors with the instrument.

If the connectors on your attenuator are straight, you can only use cables with straight connectors with the instrument.



Figure 16 Angled and Straight Connector Symbols

The above figure shows the symbols that tell you whether the optical connectors of your instrument are angled or straight. The angled contact connector symbol is colored green.

You should connect straight contact fiber end connectors with neutral sleeves to straight contact connectors, or connect angled contact fiber end connectors with green sleeves to angled contact connectors.

NOTE

You cannot connect angled non-contact fiber end connectors with orange sleeves directly to the attenuator or laser.

See "Accessories - Overview" on page 95 for further details on connector interfaces and accessories.

Making Quad-adapter Connections on the Agilent N7744/45A Multiport Power Meter




Figure 17 Connecting Optical Fibers

NOTE

Before you attach a quad-adapter or fiber, you should make sure they are clean. This is important for accurate measurements.

The fibers do not contact the optical inputs (all contact with the optical inputs should be avoided). Usually the only cleaning necessary for the optical inputs is the gentle removal of dust with pressurized, dry air.


For further information see " Maintenance and Troubleshooting - Overview " on page 97.

To use a quad-adapter to connect optical fibers to the power meter inputs:

1 Select the quad-adapter that suits the optical connector that your optical fiber is terminated with.

- **2** Attach the optical connectors to the quad-adapter. You can attach up to 4 fibers to the quad-adapter at the same time.
- **3** Make sure the number labels on the quad-adapter are at the top.
- **4** Connect the quad-adapter to the optical inputs by pressing it on until it clicks securely into place.

Removing Quad-adapter Interfaces from the Agilent N7744/5A Multiport Power Meter



Optical Output

Polarization Maintaining Fiber

A Polarization maintaining fiber (PMF) output is standard for Agilent N7711/4A tunable laser system sources. PMF is aligned to maintain the state of polarization. A well defined state of polarization helps ensure

Agilent N77xx Series

constant measurement conditions. The fiber is of Panda type, with TE mode in the slow axis in line with the connector key.




Figure 19 Polarization maintaining fiber

Electrical Connectors

There are two BNC connectors on the rear panel of your instrument. These are the Trigger Out and the Trigger In connectors.

CAUTION

The Trigger In is a TTL input. A maximum of 5 V can be applied as an external voltage to this input connector.

The Trigger Output is the only output BNC connector. This is a TTL output. Do not apply an external voltage to this connector.

LAN Interface

LAN Interface

This section explains the concept of LAN in details.

Agilent N77xx Series

Selecting a LAN Network

For the purposes of this guide, a private (isolated) LAN network is defined as a network configuration in which instrument access is a direct connection between the computer and the instrument, or to multiple instruments connected via a *dedicated* router or switch. A site (company-wide) LAN is defined as a network in which instrument access is available to many users at on-site and remote locations. The instrument's application and/or your company's IT (Information Technology) department may have guidelines that help decide the type of network (private or site) used. If a network configuration has not been determined, refer to the following considerations concerning each type.

Private LAN Considerations

Among the basic parameters of a private LAN network to consider are security, performance, reliability, and IP address availability.

Security

A private network generally involves a direct connection between the computer and the instrument, or to multiple instruments using switches or routers. Access to the instrument is limited to users on the private network, as opposed to users on a site network that could locate and access the instrument from any location.

Private networks can reduce the possibility of tests being disrupted by unplanned or unauthorized access. Code generation for test systems on a private network is often simplified, as provisions against unauthorized users may not be required.

Performance

Test systems where large amounts of data are transferred usually have faster throughput on a private network. On a site network, heavy and unpredictable LAN traffic affects each instrument (node) on the network. The impact on a test system is that repeatability is difficult to achieve as latencies are difficult to account for.

Agilent N77xx Series

Reliability

Private networks are fundamentally more reliable than site networks as they host fewer users and are less complex than site networks.

Private networks are isolated from conditions that could bring down (crash) a site network.

IP Address Availability

Every instrument (node) on a LAN (private or site) has an IP (Internet Protocol) address. Due to the expanding use of the internet, the number of site network IP addresses available is limited. By using a router with Dynamic Host Configuration Protocol (DHCP) capability on a private network, the router can assign an IP address to each instrument thus creating a sub-network (subnet) that does not consume site IP addresses.

Site LAN Considerations

For applications requiring access by many users or by users at distributed sites, a site LAN network is required. In addition to supporting multiple users, site LANs often offer the advantage of being maintained by IT departments. When using a site LAN, consult your IT department regarding LAN configuration and security issues.

Connecting the LAN Cable

LAN cables are connected between the LAN terminal on the instrument and the computer, or between the instrument and a router or switch if included in your network.

Private Network Connections

The figure below shows typical LAN cable connections for a private network.




Figure 20 Typical Private (isolated) LAN Network Connections

If the instrument is connected directly to the PC, use a LAN crossover cable, as provided with the instrument. If your computer supports Auto-MDIX or contains a LAN card with gigabit data transfer rates, the (blue) crossover cable is not required. A standard LAN cable can be used instead.

For private LAN networks that include a switch or router, use standard LAN cables for all network connections. Do not use a crossover cable.

Site Network Connections

The following figure shows a typical LAN cable connection for a site network.





Figure 21 Typical Site LAN Network Connections.

Agilent N77xx Series

On site networks, the instrument and the computer are connected directly to site LAN ports, or are connected to the site LAN through a switch. In each site network configuration, standard LAN cables are used.

The LAN LED

As the LAN connection is made, the DHCP server assigns an address and the LXI device identification proceeds you will see the following indicators:,





There may be a delay between making the LAN link (yellow status) and getting the IP address (green status). This delay may be longer if there is no DHCP server, for example when the instrument is connected directly to a PC.

The LAN Reset button

This recessed button has two functions.

• Pressing the button briefly invokes a preset of the instrument and restores default measurement settings. This is equivalent to the programming command

:SYSTem:PRESet

• Pressing and holding the button for 3 seconds will reset the LAN parameters to the factory default. This includes changing the password for modifying the configuration via LAN back to the default "agilent". This is equivalent to the following sequence of programming commands:

:SYSTem:COMMunicate:ETHernet:RESet

GPIB Interface

GPIB Interface

You can connect your GPIB interface into a star network, a linear network, or a combination star and linear network. The limitations imposed on this network are as follows:

- The total cable length cannot exceed 20 meters.
- The maximum cable length between devices is 4 meters, with an average separation of 2m between devices over the whole bus.
- No more than 15 devices may be interconnected on one bus.

Cables and Adapters

See "Cables" on page 96 for details on cables and adapters.

Connector

The following figure shows the connector and pin assignments. Connector Part Number: 1251-0293




Figure 23 GPIB Connector

Agilent N77xx Series

CAUTION	Agilent Technologies products delivered now are equipped with connectors having ISO metric-threaded lock screws and stud mounts (ISO M3.5×0.6) that are black in color. Earlier connectors may have lock screws and stud mounts with English-threaded lock screws and stud mounts (6-32 UNC) that have a shiny nickel finish.
CAUTION	It is recommended that you do not stack more than three connectors, one on top of the other. Hand-tighten the connector lock screws. Do not use a screwdriver.

GPIB Logic Levels

The instrument's GPIB lines use standard TTL logic, as follows:

True = Low = digital ground or 0 Vdc to 0.4 Vdc

False = High = open or 2.5 Vdc to 5 Vdc

All GPIB lines have LOW assertion states. High states are held at 3.0 Vdc by pull-ups within the instrument. When a line functions as an input, it requires approximately 3.2 mA to pull it low through a closure to digital ground. When a line functions as an output, it will sink up to 48 mA in the low state and approximately 0.6 mA in the high state.

NOTE

The GPIB line screens are not isolated from ground.

USB Interface

NOTE

Before connecting the instrument over USB, install the Agilent I/O Libraries Suite, which is included on a CD with the instrument.

The instrument is a USB device, with a mini-USB connector.

Powering Up the Instrument

Powering Up the Instrument

When you switch on the instrument, the LEDs on the front panel show the various stages of booting.

Standby	Power On	FPGAs loaded	Analog Board 1 loaded	Analog Board 2 loaded
Error	Yellow Error	Green Error	Green Error	Green Error
Stat	Yellow Stat	Yellow Stat	Green Stat	Green Stat
LAN	Yellow LAN	Yellow LAN	Yellow LAN	Green LAN
Yellow Pwr	Green Pwr	Green Pwr	Green Pwr	Green Pwr

Figure 24 Powering Up the Instrument

Front Panel Indicators

During operation, you may notice the following LED indicators:

	Hardware defect	
	Red	Error
LAN, USB or GPIB		
communication	Zeroing	
Green blink Stat	Yellow	Stat



Agilent N77xx Series

Connecting to the Instrument

Connecting to the Instrument - Overview

When you are first connecting the instrument, we recommend you use the Agilent Connection Expert included with the Agilent I/O libraries (available on the CD supplied with the instrument, or from www.agilent.com).

Connecting over USB

NOTE

The screenshots show the Agilent N7745A Multiport Power Meter. The same procedure applies when connecting an Agilent N7744A Multiport Power Meter, an Agilent N7751/2A Optical Attenuator and Power Meter, or an Agilent N7761/2/4/6/8A Optical Attenuator or an Agilent N7711/4A Tunable Laser System Source.

1 If it is not already running, start the Agilent Connection Expert software.





Figure 26 Agilent Connection Expert software

2 Attach the instrument to the USB. The first time you attach the instrument, follow the instructions on the screen to select the driver for this instrument.



When the instrument is connected, it is shown in the list.

Figure 27 Select the instrument in the list

3 Select the instrument in the list to see the associated tasks (in the Task Guide on the left) and information (on the right).

Finding the IP Address of an instrument

- 1 Connect to the instrument over USB, as described above.
- 2 Select the instrument in the list.

	Interface 🗙 Delete	Properties Interactive ID K Add Instrument	🖓 Refresh All 🦙 Undo
	USB Instrument - N774	Instrument I/O on this PC	ask Guide 🛛 🗙
	An instrument on the	Refresh All	asks for This Instrument
ation query was done	An identification que	E R. CND7462BNW	2 Refresh this instrument
Change Properties	VISA alias:	ST COM1 (ASRL1) COM4 (ASRL4) LAN (TCPIP0)	Change properties Send commands to this instrument
Agilent Technologies, N7745A, DEO er: Agilent Technologies	IDN string: Manufacturer:	💩 😏 USB0 177454 (USB0::2391::14104::DE00000001::0::INSTR)	Change the label
	Model code: Serial number:		🧐 alias 🎯 Ignore
0.81.xx	Firmware:		🗙 Delete
USB0::2391::14104::DE00000001::0 usb0[2391::14104::DE00000001::0	VISA address: SICL address:		General Tasks
	Address check:		2 Refresh al
Yes Idor Agilent Technologies	Auto-identify:		Nore Information
Idor: Agrene reciriologies	USB driver vendor:		iore anormadori
			How do I get drivers?
			Where can I find programming samples?
1		<	Where can I find programming samples?

Figure 28 Select the instrument

Agilent N77xx Series

3 Click on "Send commands to this instrument".

The default command is the *IDN? query.

🕌 Agilent	Interactive I	D - CONNE	CTED TO USE	30::2391::14	¥1 🔳	
Connect In	teract Help					
(2) Stop	Device Clear	101 Read STB	SYST:ERR?	<u>ラ</u> Clear History	📰 Options	
Command:	"IDN?			~	Commar	nds 🕨
, i	Send Command ession History:	Read Re	sponse Sei	nd & Read		
* Conne	cted to: USI	30::2391:	:14104::DEO	0000001::0:	: INSTR	< >

Figure 29 *IDN? query

4 Enter the command

:SYSTem:COMMunicate:ETHernet:IPADdress:CURRent?

then click on [Send & Read].

The instrument returns its current IP address.

🐴 Agilent	Interactive IO - CONNECTED TO USB0::2391::141 🔳 🗖 🔀
Connect In	teract Help
8	p 🔟 🖳 词 📰
Stop	Device Clear Read STB SYST:ERR? Clear History Options
Command:	SYSTem:COMMunicate:ETHernet:IPADdress:CURRent? 🔽 Commands 🕨
(Send Command Read Response Send & Read
Instrument S	ession History:
-> :SYST	cted to: USB0::2391::14104::DE00000001::0::INSTR American Ame American American Am American American Americ American American Ame

Figure 30 :SYSTem:COMMunicate:ETHernet:IPADdress:CURRent?

Connecting over LAN

NOTE

The screenshots show the Agilent N7745A Multiport Power Meter. The same procedure applies when connecting to other N77xx instruments like attenuators and laser sources.

- 1 Make sure the instrument is connected to the LAN, and that the LAN LED on the front panel is green. See, "LAN Interface" on page 28 for more information to connecting over LAN.
- **2** If it is not already running, start the Agilent Connection Expert software.
- **3** Select "LAN (TCPIP0)" in the list.




Figure 31

- 4 Click on "Add an instrument" in the Task Guide.
- **5** Wait for the scan to finish.
- **6** Select the instrument in the list.

Agilent N77xx Series




Figure 32

7 To show the description of the N775xA or N776xA, check the box for *IDN? query and click on "Find Again".

Auto Discover li	Find Add Address Manually add a known Padress or hostname	Search Explore a network address range network	aments	
Select	Address (IP, MAC, Host)	Description	Instrument Web Page	Automatically find and identify local instruments.
	134.40.93.152 00-30-d3-09-f9-de a-n4385a-50pp06.germany.agilent.com	(No description available)	Web Page	A local instrument is one on the same subnet as any of the network interfaces in
	134.40.93.243 00-30-d3-0e-5f-49 b-n4387a-50rd01.germany.agilent.com	(No description available)	Web Page	the computer.
	134.40.93.247 00-30-d3-0e-be-84 loxi2.germany.agilent.com	P1K Example (000000000)	Web Page	
	134.40.95.120 00-30-d3-09-59-44 dhcp-bbn6-1.germany.agilent.com	P1K Example (0000000000)	Web Page	
	134.40.95.141 00-30-d3-0e-7f-b5 embed-ppmc5.germany.agilent.com	P1K Example (0000000000)	Web Page	Search on this page for:
V	134,40,95.155 00-30-d3-0e-7f-b4 dhcp-bbn6-7.germany.agilent.com	P1000	Web Page	Allow "IDN? query
				Find Again

Figure 33

- 8 Select the instrument in the list to see the associated tasks (in the Task Guide on the left) and information (on the right).
- For the N774x Multiport power meters, click on the "Instrument Web Interface" button to control the instrument directly.



Figure 34

If it is not already installed, you will need to install a Java run-time engine to use this page.

Basic configuration and functionality are available from this web page.

- For the attenuators, and as an alternative for the power meters, you can control the instrument from the Agilent N77xx Viewer.
 - Connecting using the N77xx Viewer is described in .
 - Making power measurements with the N77xx Viewer is described in "Controlling the Optical Multiport Power Meter" on page 49.
 - Controlling attenuation with the N77xx Viewer is described in .

More advanced functionality is available through programming with SCPI commands (described in the Programming Guide available on the CD supplied with the instrument or from www.agilent.com), or using the 816x VXI Plug&Play driver, version 4.3 or higher. For the N774x Multiport power meters you can also use the IVI-COM driver.

Agilent N77xx Series

Claims and Repackaging

Claims and Repackaging

If physical damage is evident or if the instrument does not meet specification when received, notify the carrier and the nearest Agilent Technologies Sales/Service Office. The Agilent Technologies Sales/ Service Office will arrange for repair or replacement of the unit without waiting for settlement of the claim against the carrier.

Return Shipments to Agilent Technologies

If the instrument is to be shipped to an Agilent Technologies Sales/ Service Office, attach a tag showing owner, return address, model number and full serial number and the type of service required.

The original shipping carton and packing material may be reusable, but the Agilent Technologies Sales/Service Office will provide information and recommendations on materials to be used if the original packing is no longer available or reusable. General instructions for repackaging are as follows:

- Wrap instrument in heavy paper or plastic.
- Use strong shipping container. A double wall carton made of 350pound test material is adequate.
- Use enough shock absorbing material (3 to 4 inch layer) around all sides of the instrument to provide a firm cushion and prevent movement inside container. Protect control panel with cardboard.
- Seal shipping container securely.
- Mark shipping container FRAGILE to encourage careful handling.
- In any correspondence, refer to instrument by model number and serial number.

Deleting user data

If you need to delete all your logged data and user configurations, that is to reset the instrument completely:

- **1** Press and hold the Reset/standby button until the instrument starts to reboot.
- 2 When the instrument reboots, press the Reset/standby button once more briefly, to confirm you want to delete all the stored data.

N77xx Viewer Common User Interface Functions

Connecting to an Instrument

When you first launch the N77xx Viewer, it is not connected to any of your instruments.

- 1 Click on the "Connect" icon at the top right of the window.
- 2 Select the instrument to which you want to connect from the list.

Agilent N77xx Viewer - C	onnect	
Select an instrument from the TCPIP0::134.40.95.213::inst0 TCPIP0:A-N7764A-300075.c USB0::2391::14104::DE4920 USB0::2391::14104::DE0000	t:INSTR termany.agilent.com::inst0::INSTR 0042::0::INSTR	
Scan Again	Place a shortcut on the Desktop	OK Cancel
Instrument Data		
VISA Address	TCPIP0::A-N7764A-300075.germa	ny.agilent.com::inst0::INSTR
IP Address	134.40.93.232	
MAC Address	00-30-D3-11-D2-75	
Identification (*IDN?)	Agilent Technologies, N7764A, DE	49300075, 0.99
Manufacturer	Agilent Technologies	
Model Code	N7764A	
Serial Number	DE49300075	Times of the set
Firmware	0.99	Lange Article
		Configure Interfaces

Figure 35 Connecting to an Instrument

Agilent N77xx Series

3 The list shows instruments that have been given a VISA address.

If your instrument is not in the list,

- 1 Use the Agilent Connection Expert to check your instrument has a VISA address
- 2 Click on the "Scan Again" button.

Once communication has been established to the instrument, its details are shown.

4 If you need to change the connection configuration for an instrument on the LAN, click on the "Configure Interfaces" button.

You can only configure the LAN parameters before the instrument is connected.

Interface Param	eters	
	Auto IP	
IP Address	manual 127.0.0.1	current 134.40.93.232
Subnet Mask	255.255.255.0	255.255.248.0
Default Gateway	0.0.0.0	134.40.88.1
DNS Server	-	
Host Name	A-N7764A-300075	1
Domain Name		
Net BIOS		
Ethemet Conne	iction Monitoring	
Description		
TCP Keep Aller		
TCP Keep Alive Tin	ne 🛛	
GP18		
SICL Interface		
GPIB Address	20	
Security		
Password		
	0	K Cancel

Figure 36 Interface Parameters

- Make sure both DHCP and Auto IP are not selected.
- You can now edit the IP parameters.
- **5** For instruments you use regularly, you can add a link on your desktop by putting a check mark in the box to "Place a shortcut on the Desktop".
- **6** Click on the "OK" button to finish the connection.

Viewing Measurements

There are two ways you can view results for a power meter or attenuator.

Viewing the measurement on a single channel

1 Click on the "New" button to open a window with the current measurement and parameter settings.

Attenuator 1	->2 🔳 🔀	Laser	
-10 Output	.002 dBm		9000 THz
C set P set Power Control Wavelength Avg Time C Offset P Offset Trigger In Trigger Out	4.038 dB -10.000 dBm On 1547.5 nm 2 ms 0.000 dB 0.000 dB Ignore Disabled	Output Wavelength Channel f ₀ Δf Ch. Spacing FM FM Deviation AM Level	0.000 mVV 1554.134 nm -2 193.1 THz 0 GHz 100 GHz Off 0 GHz
7		,	

Figure 37 Viewing the measurement on a single channel attenuator/power meter (left) or tunable laser (right)

Viewing all measurements

1 Click on Overview to open a window with the measurements on all channels.

owermeter 1		Powermeter 2		Powermeter 3		Powermeter 4	
1	.767 pW	1	.777 pW	1	.374 pW	1	.521 pW
√avelength	1550 nm	Wavelength	1550 nm	Wavelength	1550 nm	Wavelength	1550 nm
∖vg Time	200 ms	Avg Time	200 ms	Avg Time	200 ms	Avg Time	200 ms
Range 💼	-30 dBm	Range 😐	-30 dBm	Range 💼	-30 dBm	Range 💶	-30 dBm
Ref. Mode	Off						
Ref. Channel	0.000 dB						
lef. Value	1.000 uW	Ref. Value	1.000 uW	Ref. Value	1.000 uW	Ref. Value	1.000 uW
alibration	0.000 dB	Calibration	0.000 dB	Calibration	0.000 dB	Calibration	0.000 dB
rigger In	Ignore	Trigger In	Ignore	Trigger In	Ignore	Trigger In	Ignore
rigger Out	Disabled	Trigger Out	Disabled	Trigger Out	Disabled	Trigger Out	Disabled

Figure 38 Viewing all measurements

Saving and Recalling Configurations

Saving a configuration

You can save the current measurement configuration to either of two settings stored on the instrument.

- **1** In the File menu, select Save.
- 2 Enter or select the setting to which you want to save it.

Select the	e Setting number
Minimur	n: 1; Maximum: 2
Ħ	÷
	OK Cancel

Figure 39 Saving a configuration

3 Click on the "OK" button.

Recalling a configuration

- 1 In the File menu, select Load.
- 2 Enter or select the setting you want to retrieve.
- **3** Click on the "OK" button.

Controlling the Instrument Refresh

1 In the File menu, select Update.

In the dialog, you can :

Instrument Update	Cycle Settings	
Update enab	led	
Readings	350	ms
Polled Values	2000	ms
Default	OK Can	cel

Figure 40 Controlling the Instrument Refresh

- Stop and start the display of all changes. Make sure to put a check in the box for "Update enabled", to enable the N77xx Viewer to display the status and results from the instrument. If there is no check in the box for "Update enabled", the display is not updated with the replies to queries sent to the instrument.
- Set the update interval for "Readings" to control how often the result is read from the instrument and the display updated.
- Set the update interval for "Polled Values" to control how often the measurement parameters are read from the instrument and the display updated. This is useful if, for example, more than one person can control the instrument.
- **2** When you have set all the values, click on the "OK" button.

Agilent N77xx Series

1 Getting Started



Power Measurement

This chapter describes how to measure optical power using the Agilent N7744A/N7745A Optical Multiport Power Meter or the power meter channels of the Agilent N7751A/N7752A Optical Attenuator and Power Meter.

- "Controlling the Optical Multiport Power Meter" on page 49
- "The Power Value" on page 53

Controlling the Optical Multiport Power Meter

There are different possibilities to control the instrument:


- LAN Connection: [Instrument Web Interface], for the Agilent N7744A/N7745A Optical Multiport Power Meter only.
- LAN / USB / GPIB Connection: [N77xx Viewer]
- LAN / USB / GPIB: Programming with SCPI commands (described in the Programming Guide available on the CD supplied with the instrument or from www.agilent.com), or using the 816x VXI Plug&Play driver. For the Agilent N7744A/N7745A Optical Multiport Power Meter, you can also use the associated IVI-COM drivers.

To connect to the instrument see " Input and Output Connectors -Overview" on page 21.

To display the Web Interface of the Multiport Power Meter please either

- Press the "Instrument Web Interface" Button within the Connection Expert, or
- Enter the instrument's IP Address or DNS Address in the Address Line of your browser.





The browser starts with the screen shown in the figure, below.

Figure 41 The Welcome Page

Select the "Browser Web Control" tab on the left.

P Observe Only Whene Full Condition	Carpine Preser		lingby Nin Mar Correr	with:	LXI	
In Plantag (former)		×	1 Monte			
	0.812	mW		SH	Segn	
1 2	3.717	pW			Port #	No 1 💌
# 3	1.730		Paran	vetere	Ri	atonasce
	1.526	PW		u		AI.
# 4		Wg	Wavelength (http://		Reference Mode	C Of
# 5	2.995	₽W	Retaging Time 1	0.7 9 😽	Reference Channel	0.000 08
# 6	2.864	pW	Provec Unit 1	W/. 🖌	Reference Yalue	1 000 µ//
. 7	1.994	DM	Power Rainge 🗍			Displayin fi
# 8	3.891	₽₩	Range Moda (automatic 💌	Catanatoo	17 0.000 at
Pause		Close	Moste	laneous		riMan
-				N.		61
Billiping Monthéon		*	Power Sigits	3	MinMar Blaze	r on w
: AP -	Misinam / Misinum	Current			NiniMar Może	T Wedaw
COLUMN STORE	1 100 W 1000 W	TO SHOW THE	[#1 □ Qued 1	F Qued 2	Data Points	F 1060
1 1 1 100 K 1 7 1 100 K	8 800 W K.000 W 8 800 W K.000 W	8,000 H 1,000 H	Lero			Fesel
1 1 1.150 T	1 000 T 1 000 T	8 000 T 8 000 T 1 000 H	and the second	i manual i		
A 1 4 4 5 4	8 000 W 8 000 W	1.000 W	Solu 1 - Porol	Freed		

Figure 42 The Browser Web Control

How to Change the Settings of the Optical Multiport Power Meter

To change the settings of the Optical Multiport Power Meter from the Browser Web control window, you must select the "Allow Full Control" button. Now you are able to change the settings (see the figure, below) of the Optical Multiport Power Meter. The changes are made when you press the "Apply" button. If you selected the "All" checkbox, the setting is applied to all ports.

To quickly and reproducibly set the measurement parameters, use the Measurement Settings buffers. These let you save and recall sets of settings with single programming commands. For further details, refer to "Handling Measurement Settings -

The :CONFigure:MEASurement:SETTing subtree" in the Programming Guide.

		Setti	nga				
			Port#		No.1	~	
P	arameter	8	Re	ference			
	68			AL			
(Yavelength (nm)	Г	1556.0	Reference Mode	T.	or		v
weraging Time	E	0.2 8 😽	Reference Channel	E.	0.080	dÐ	
Power Unit	F	W 💌	Reference Value	F	1.000	10 W	٣
Power Range	F			Г	Displa	ta Ref	etence
Range Mode	Г	automatic 💌	Calibration	Г	0 000	dÐ	
Power Digits	AR F	a 💌	Min/Max Slate	All F	or	×	
			Um/Mai: Mode	F	Windo	et :	*
	d1 [Quad 2	Data Points	Г	1000		
Zi	era			-	R	eset	
Save 1 💌 R	ecali	Preset					
	los more de	in all	OK	Close	1	ophy.	- 11

Figure 43 Settings of the Optical Multiport Power Meter

Agilent N77xx Series

The most convenient way to control one or more Multiport Power Meters is using the N77xx Viewer, which is a control and display graphical user interface.

😹 N7744A DE49200	042 - N77xx Viewer	
Instrument View H	telp	
the second residence of the second	efresh Overview New MinMax Help	
Powermeter 1 Powerr	neter 2 Powermeter 3 Powermeter 4	
[]		
	-260.271 pW	To Ref
Settings		
Wavelength	1550 nm 🔣	
Avg Time	200.000 ms 💌 🔣	
Power Unit	Watt 💌 🖂	
Power Range	-30 dBm 👻	
Range Mode	Auto 💌 🚭	
Reference and	Calibration	
Ref. Mode	Off 🔽	
Ref. Channel	0.000 dB 🔣	
Ref. Value	1.000 W 💙 📩	
Calibration	0.000 dB	
Display Digits	3	
	nt Technologies, N7744A, DE49200042, 1.12 STR Agilent N77xx Viewer - Version 1.2.0.0	

Figure 44 The N77xx Viewer

The Viewer is available to download from the product webpage which you can access at

http://www.agilent.com/find/mppm

Then select the tab "Technical Support" followed by "Drivers & Software".

earch: Drivers & Software	~	Go		∰ Store
show all images #	Product Stat	upport: tical Multiport Power Meter tus: Currently Orderable Currently Supported les: Contact Agilent for details	8 <u>Need</u>	Support Now [+] Feedbac
Drivers & Software Refine the List	rivers & Software	Application Notes Specifications FAQs more		t Overview
By Class		Title/Description	Date	Туре
VX (2) VXIplug&play (2) By Application @ Digital & Analog Design	n (1)	THEP/DESCRIPTION <u>816x VXI Plug&Play Driver</u> This popular driver greatly simplifies programming applications with the 816x Lightwave Measurement System, especially for coordinating tunable lasers with optical power meters for swept spectral measurements.		Driver
By Type of Content				

Figure 45 Technical Support

On this webpage you will also find additional software tools like:

- Firmware update utility
- Firmware update file
- VXI Plug&Play Drivers
- IVI Drivers

How to Measure Power

The Power Value

For each port, the optical power is measured and returned or displayed. The power reading represents the average optical power over the selected averaging time. The power meter can continuously self-trigger

Agilent N77xx Series

new measurements or be configured to wait for external software or hardware triggers.

How to Set the Number of Digits

When the power is displayed with the Web Interface or the Viewer program, you can set the maximum number of digits that are used in power measurement. This is the maximum number of digits after the decimal point.

How to Set the Power Unit

The Power Unit setting allows you to select either W, dBm or dB as the units in which power is displayed or returned. For returning large amounts of logged data, the fastest data transfer is achieved using W units.

What are the Power Units?

Watts (W) are the SI unit for power measurement.

You can also measure power in dB or dBm. Values returned in these units are derived from measurement in Watts.

By selecting dBm, the following calculation is made:

$$P_{dBm} = 10\log \frac{P_{input}(W)}{1 \times 10^{-3}(W)}$$

Where,

 \boldsymbol{P}_{dBm} is the power value returned in dBm, and

P_{input} is the input signal level in Watts.

Power, in units of dBm, is measured relative to 1 mW. It is an absolute power measurement.

Power, in units of dB, is measured relative to a particular reference power value. For information on selecting this reference value, see How to Input a Reference Level.

The following calculation is made:

$$P_{dB} = 10\log \frac{P_{input}(W)}{P_{ref}(W)}$$

Where,

 P_{dB} is the power value returned in dB,

P_{input} is the input signal level in Watts, and

P_{ref} is the chosen reference power value in Watts.

How to Set the Calibration Offset

This is a calibration offset that you can enter to compensate for external optical circuitry, for example, an attenuator, switch, or a 3 dB coupler. The calibration offset, [Cal], is automatically subtracted from the input signal, as expressed in dB.

$$P_{measured}(dBm) = P_{input}(dBm) - CAL(dB)$$

where,

P_{measured} is the adjusted value of the signal read in dBm,

P_{input} is the input signal level in dBm, and

CAL is the calibration offset in decibels (dB).

How to Set the Reference Level

How to Input a Reference Level

dB results are shown relative to a reference level. Setting, or changing, the reference only affects results that are returned in dB.

You can set the reference level in units of dBm or Watts. Depending on which is selected the following equations are used to calculate the power level in dB:

 $P_{dB} = P_{measured}(dBm) - REF(dBm)$

or

$$P_{dB} = 10\log \frac{P_{measured}(W)}{REF(W)}$$

where,

 P_{dB} is the returned relative power,

 P_{measured} is the absolute power level (see "How to Set the Calibration Offset" on page 55), and

Agilent N77xx Series

2 Power Measurement

REF is the reference level.

How to Set the Reference Value to the Current Power Value

In addition to entering a new reference value, you can change the reference value to the current power value from the N77xx Viewer or by programming.

This takes the input power level and stores it as the reference. Setting the reference only affects results returned in dB.

The power value in dBm or Watts is stored as the reference, that is:

 $REF = P_{measured}$

where,

REF is the reference, and

 P_{measured} is the absolute power level (see "How to Set the Calibration Offset" on page 55).

How to Reference Another Power Measurement Port

You can use the reference mode to reference to another port's current power level. If you choose this mode, the reference of the port continuously updates to the current power value measured by the other port.

The following equations describe how the power in dB is calculated if you reference another port:

 $P_{dB} = P_{measuredA}(dBm) - P_{measuredB}(dBm) - REF(dB)$

or

$$P_{dB} = 10\log \frac{P_{measuredA}(W)}{P_{measuredB}(W)} - REF(dB)$$

where:

 P_{dB} is the relative power,

 $P_{measuredA}$ is the absolute power level (see "How to Set the Calibration Offset" on page 55) measured by the current port, and

 $P_{measuredB}$ is the absolute power level (see "How to Set the Calibration Offset" on page 55) measured by the port you choose to reference, and *REF* is the reference level in decibels (dB).

How to Set the Wavelength

This is the wavelength value of the signal to be measured. The responsivity of the Power Meter varies with wavelength.

For accurate absolute power measurement, you need to input the wavelength of the optical input.

For relative power measurements, this setting is not critical.

How to Remove Electrical Offsets (Zero)

Optical Power Meters measure optical power by converting optical power to electrical current, and then measuring electrical current. An electrical offset is electrical current that is always present, even if there is no optical power input. If electrical offsets are not removed, they affect the accuracy of power measurement, especially at low optical power levels.

Performing a zero sets the zero power level to the average electrical offset level for the current environmental conditions.

NOTE	The environmental conditions and the temperature of the instrument NOTE affect
	electrical offset.

For the best results you must:

- Allow the instrument time to acclimatize (around 24 hours).
- Allow the instrument time to warm up (around 20 minutes).
- Make sure that the optical input is not receiving any light or use the zeroing adapter. If you are using multi-mode fiberoptic cable, you must disconnect the cable and cover the input to the Power Meter to perform a zero.

It is good practice to perform a zero before making any sensitive measurements.

Logging data over a wide dynamic range is a typical case that requires a valid zero setting.

1 In the File menu, select Zeroing.

Agilent N77xx Series

2 Power Measurement




Figure 46

2 Select the channel you want to zero from the drop down list and click on the "Port" button.

OR

Click on the "All" button.

While the zero is running, the traffic light flashes amber.

If you receive a zero error, the zeroing operation has failed because the Power Meter received light. The most common reason for zeroing to fail is if:

- a source is connected to the Power Meter's input connector,
- the fiber connected to the Power Meter's input connector is collecting light, or
- the Power Meter receives ambient light because the input connector is uncovered.

When zeroing is finished, the traffic light is green. You can now proceed with your measurements.

When environmental conditions are stable, the zero setting will remain valid over long times and the same settings will be active after rebooting the instrument.

How to Choose the Range Mode

You can choose either of two ranging modes:

- *<Auto>*, the auto-ranging mode, ensures that the result has a returned value between 9% and 100% of full scale. The default state is for automatic ranging to be enabled.
- *<Manual>*, which allows you to set a user-defined range.

How to Set the Range

The range setting determines the power levels that can be measured. The setting should be high enough for the applied power, but lower settings will provide a lower noise level for weak signals.

If you choose auto-ranging mode, this parameter can not be set. The Range parameter returns the automatically-determined range.

If the measured power is more than 3dB greater than the range setting, the power cannot be measured and an overflow is indicated by the returned power value. This means that the measured power is greater than the Upper Power Limit. You must increase the [*Range*] value.

More more details, see the following table.

Upper Power Limit

The following table shows the upper power limit at the available power ranges. The lowest measurable power in each range is determined by the noise level, as indicated in the specifications.

Table 4	Upper Power Limits for Various Power Ranges			
[RANGE]	Upper Power Limit			
10 dBm	10.000 mW	10.000 dBm		
0 dBm	1.9999 mW	3.000 dBm		
-10 dBm	199.99 µW	-7.000 dBm		
-20 dBm	19.999 µW	-17.000 dBm		
-30 dBm	1.9999 µW	-27.000 dBm		

Agilent N77xx Series

How to Set the Averaging Time

This is the length of time over which a signal is averaged to produce a measured value. Longer averaging times increase the accuracy and improve the noise rejection.

The averaging time of the N7744A and N7745A can be set over a wide range from 1 μ s to 10s. Up to 2 ms, the setting can be varied with a typical resolution of 2 μ s, and for higher averaging times the setting resolution is about 0.1%. This flexible setting time supports optimal synchronization with the time scale of the signals.

The high sampling rate with short averaging times is supported by the high cutoff frequency, as shown in the specifications. These bandwidths are reduced in steps, when longer averaging times are used to improve the noise performance while maintaining signal fidelity.

How to Set Automatic Gain for Averaging Modulated Signals

A fast Auto Gain switching is used to achieve the high dynamic range. In some applications, like determining the average power of signals modulated with periods shorter than the averaging time, this switching produces irregular readings. In such cases, the Auto Gain should be turned off, as described in the programming guide. This setting takes effect after setting the averaging time on the port directly or via function parameters.

How to Choose the MinMax Mode

MinMax mode measures the incoming power and returns the minimum value measured, and the maximum value measured. The difference between these values, delta P, is returned in place of P, the power value. This mode is intended principally for polarization dependent measurements, but can be used for other types of measurement.

Powermeter 1		Powermeter 2		Powermeter 3		Powermeter 4		
ΔΡ	1.052 dB	ΔΡ	1.116 dB	ΔΡ	2.088 dB	ΔΡ	1.502 dB	
Pmin	-87.754 dBm	Pmin	-87.252 dBm	Pmin	-89.833 dBm	Pmin	-89.355 dBm	
Pmax	-86.702 dBm	Pmax	-86.135 dBm	Pmax	-87.745 dBm	Pmax	-87.853 dBm	
MinMax Mode	Window	MinMax Mode	Window	MinMax Mode	Window 🔽	MinMax Mode	Window 👻	
Data Points	1000							
Wavelength	1550 nm							
Avg Time	200 ms							
Range	-30 dBm							
Ref. Mode	Off							
Ref. Channel	0.000 dB							
Ref. Value	-30.000 dBm							
Calibration	0.000 dB							
Stop		Stop		Stop		Stop		
<u> </u>				<u> </u>		<u> </u>		

Figure 47 How to Choose the MinMax Mode

The Min/Max mode can be activated from both the Web Control Interface and the N77xx Viewer, as well as with direct programming described in the Programming Guide.You can choose two modes of operation from the MinMax mode menu:

• *<Continuous>* mode, which compares each new measured value with the maximum and minimum values so far, and replaces them as necessary.

This mode is useful for measuring the Polarization Dependent Loss (PDL) of a component. Run the application while sweeping the polarization of the source applied to the component.

• *<Refresh>* mode, which records the minimum and maximum values until N samples have been compared. After N samples the minimum and maximum values are reset and a new MinMax measurement begins.

You can use Refresh mode, for example, when you are searching for or setting the position of minimum PDL.

Refresh Mode	▲ N Samples ↓ Time for display to update	The length of the lines represents the size of the buffer at the time of update
	·	
	_	

Figure 48 Refresh Mode

Agilent N77xx Series

2 Power Measurement

N O T E	Window mode is not currently implemented in the N774xA and N775xA. For compatibility with existing software it can be selected but functions the same as Refresh mode.
NOTE	The minimum sampling period in Min/Max mode, dependent on the rate of commands and number of active ports, is about 2 ms.

How to Use External Triggers

The Agilent N7744A / N7745A Optical Multiport Power Meter lets you trigger the instrument to perform tasks and to output trigger signals to external measurement instruments using the rear panel BNC connectors.

1 In the File menu, select Trigger.

	e Trigger Configuration	Default	~					Res	et All
Configuration	n Choices								
🥑 <dis< th=""><th>abled>: disables back</th><th>plane BNC trigge</th><th>r connectors.</th><th></th><th></th><th></th><th></th><th></th><th></th></dis<>	abled>: disables back	plane BNC trigge	r connectors.						
<del< th=""><th>fault>: input and outpu</th><th>t trigger connect</th><th>ors are enabled to pa</th><th>ss trigger signals.</th><th></th><th></th><th></th><th></th><th></th></del<>	fault>: input and outpu	t trigger connect	ors are enabled to pa	ss trigger signals.					
<pa:< td=""><td>ss Through>: input cor</td><td>nnector is enable</td><td>d to pass trigger signa</td><td>ls to instrument a</td><td>and to output c</td><td>onnecto</td><td></td><td></td><td></td></pa:<>	ss Through>: input cor	nnector is enable	d to pass trigger signa	ls to instrument a	and to output c	onnecto			
	2 3								
nannei i rig <u>c</u> Powermeter	ger Functions		~	Powermeter	2				
		0.0000000000000000000000000000000000000	Powermeter 2			1280	Powermeter 4		
Trigger In	Single	Trigger In	Single 🔽	Trigger In	Complete	~	Trigger In	Ignore	
Trigger Out	Meas. Started 💉	Trigger Out	Avg Time Over ⊻	Trigger Out	Disabled	~	Trigger Out	Disabled	1
<cor< th=""><th>ore>: input triggers are mplete>: an input trigge gle>: an input trigger w</th><th>er will trigger a co</th><th>mplete measurement</th><th>function to be pe</th><th></th><th></th><th></th><th></th><th></th></cor<>	ore>: input triggers are mplete>: an input trigge gle>: an input trigger w	er will trigger a co	mplete measurement	function to be pe					
Trigger Out (Choices								
	abled>: the output trig	ger mode is disab	oled.						
💡 <dis< td=""><td></td><td></td><td></td><td>15.</td><td></td><td></td><td></td><td></td><td></td></dis<>				15.					
- 	as. Started>: a trigger								

Figure 49 Trigger Configuration

2 Set the trigger configuration

By default, the trigger level at the external trigger connector is active high. This means when a trigger rises above the high TTL level (rising edge), a trigger is accepted.

You can select three modes of triggering from the trigger configuration menu:

- <Disabled>, if you do want to disable triggering.
- <Default> to enable the trigger connectors.
- <Pass Through> if you want an input trigger to automatically generate an output trigger. This lets you trigger another instrument almost simultaneously.
- **3** Set the input trigger functions

The power meter ports can be configured individually to respond to an input trigger by:

- <Ignore> The trigger is ignored.
- <Single Measurement> The trigger causes a single power measurement to be made.
- <Complete Measurement> The trigger will cause an active measurement function like logging to begin and continue until completed. Please note, logging and stability functions may only be available on the Agilent N7744/5A Multiport Power Meter.
- **4** Set the output trigger functions

The power meter ports can be configured individually to generate output triggers for the following events:

- <Disabled> No triggers are generated.
- <Meas. started> An output trigger is generated when a power measurement begins.
- <Avg Time Over> An inverted trigger is generated for each measurement sample.

For further information on triggering, please refer to the Agilent N77-Series Programming GuideGuide, in the chapter "Triggering - The TRIGger Subsystem".

Agilent N77xx Series

2 Power Measurement

www.valuetronics.com

Agilent N77xx Series



Power Meter Applications

This chapter explains how to set up and perform the following applications using the N7744A or N7745A Optical Multiport Power Meter.

The Agilent N775xA power meters do not support the Logging and Stability Applications.

- The Logging application logs power measurements. The Logging application performs measurements consecutively without the possibility of any pause between measurements.
- The Stability application logs power measurements. The Stability application differs from the Logging application because:
 - the period time, the time it takes from the start of one measurement period until another measurement period starts, may be greater than or equal to the averaging time.

These applications can be accessed with programming commands as described in the Programming Guide.

- "The Logging Application" on page 65
- "The Stability Application" on page 67

The Logging Application

The Logging Application logs a series of power measurements for a number of power meter ports.



3 Power Meter Applications




Figure 50 Example Logging Application

The most important parameters for the logging application are:

- *[Samples]*, this is the number of samples that will be performed during the logging application, and
- *[AvgTime]*, the averaging time is the length of time over which a signal is averaged. Longer averaging times increase the accuracy and improve the noise rejection. Longer averaging times also improve sensitivity.

The total time for a logging application is the number of samples multiplied by the averaging time. The power range is constant during a logging application and autoranging is disabled.

The Logging application differs from the "The Stability Application" (see, "The Stability Application" on page 67) because all power measurements are performed consecutively without the possibility of any pause between measurements, as shown in the figure, above.

While performing the logging application all used ports of the N7744A / N7745A Multiport Power Meter will have the same averaging time and number of samples. Setup parameters for each of the ports should be set before the logging is started on any port.

The logging function is often used together with singlemeasurement input triggering for synchronized measurement with another instrument, especially for spectral measurements with a tunable laser.

Continuous Logging Using Data Buffers

Enhanced logging functionality is available using two data buffers for each port.

The MPPM ensures there is no time lag between logging and data availability by using a data buffer with a capacity of 1M samples for each channel. You can read the most recent results out of one buffer out while the next logging measurement is filling the other buffer.

You can use this to speed up applications with repeated logging operations. It is especially valuable for monitoring signals over extended periods to detect transient events.

To use the data buffers, set the logging function to perform a fixed or indefinite number of "LOOP"s, as described in the Programming Guide.

- When LOOP = 1, default logging behavior is selected, that is measurements are made and the results are written to the first buffer.
- If LOOP = 0, logging continues, writing the results into alternate buffers, until you stop it.
- When LOOP = n>1, then upon completion of a logging measurement, another is started and the results are written in rotation to the two buffers (first to buffer A, then buffer B, then buffer A, and so on).

This is repeated until n logging measurements have been performed. For example, if n=2, both buffers are filled, and you can upload a total of 2M samples for the channel.

You can also use input triggers to begin individual loops or to begin multiple loops with an initial trigger.

The "Index" value is updated to indicate that data from each LOOP is available for readout.

The Stability Application

The Stability Application logs a series of power measurements for a number of power meter ports.




Figure 51 Example Stability Application

The Stability application differs from the "The Logging Application" (see,) because you set the averaging time, the total time, and the period time after which each new measurement is started. Period time must be greater than or equal to the averaging time you set.

This leaves the following two possibilities:

- The averaging time is less than the period time, there is a dwell time between the commencement of each power measurement, as shown in the figure, above.
- The averaging time is equal to the period time, all power measuremant are performed without any dwell time.

ΝΟΤΕ

The total time concludes after the last dwell time period, not after the last averaging time period.

The Stability application enables you to make periodic power measurements over a long time (277 hours (11 days and 13 hours)).



Definition of Terms

(Applicable) fiber type


Connection type for which the specifications and characteristics apply (if not differently stated).

Absolute wavelength (frequency) accuracy

The maximum difference between the displayed wavelength (frequency) and the actual wavelength (frequency) of the tunable laser source. Wavelength is defined as wavelength in vacuum.

Attenuation

Difference (in dB) between "Total loss" (with any attenuation setting, see "Total loss" on page 86) and "Insertion loss" (i.e. with attenuation set to zero, see "Insertion loss" on page 73) of the attenuator.



Attenuation setting mode

Operating mode where the user sets the desired attenuation. The power control function is deactivated in this mode (see also "Power setting mode" on page 80).

Attenuation accuracy (uncertainty)

The maximum possible difference (in dB) between the displayed "Attenuation" (see, "Attenuation" on page 69) and actual "Attenuation" (see, "Attenuation" on page 69).




Figure 53 Attenuation accuracy (uncertainty)

Conditions: Attenuator set to the wavelength of the source. Input polarization, temperature range, maximum power and input mode (multimode attenuators only) as specified. Other conditions as specified.

Measurement: With Fabry-Pérot laser source of applicable wavelength, polarization scrambler (if applicable), and optical power meter.

When using a polarized source, the attenuation uncertainty increases by \pm half the specified polarization dependent loss of the attenuator (if specified).

Averaging time

Specifies the period during which the power meter or monitor power meter takes readings from the photodetector. The displayed power result is the (arithmetical) average of these readings.

NOTE

Constant operating conditions

This generally includes constant values of temperature, humidity, wavelength, input power level, polarization state and mode distribution, if the quantity is not explicitly subject to variation.

Data logging capability

Maximum number of power readings per port that can be captured during a logging operation.

Drift (dark)

Specifies the maximum deviation from zero of the measured power with no optical input (dark) over time, excluding "Noise pp (dark)" on page 74, power expressed in Watts.

Measurement: From the dark measurement result P_i , the 3rd order polynomial fit F_i is calculated. The fit incorporates the least square-sum difference. The Drift D is calculated by

 $D = max_i \{F_i\}$ or $mini \{F_i\}$, whichever is worse.




Figure 54 Drift (dark)

Conditions: All ports dark. Observation time as specified. Other conditions as specified.

Dynamic range (logging mode)

For a fixed power meter range, dynamic range D specifies the difference between the maximum displayable value P_{max} in dBm in the corresponding power meter range and the minimum measurable power P_{min} in dBm.

Conditions: As specified.

Measurement: P_{max} is 3 dB above the nominal power meter range value (that is, in the +0 dBm power meter range P_{max} is +3 dBm), except for the +10 dBm range with P_{max} = +10 dBm. P_{min} , dBm is calculated as the root-sum-square (RSS) of the specified "Noise 2 sigma (dark)" (see, "Noise 2 sigma (dark)" on page 75 and "Drift (dark)" (see, "Drift (dark)" on page 71. RSS is calculated in linear units (Watts).

External analog modulation - modulation depth

Specifies half the peak-to-peak optical power change divided by the average optical power for a sinusoidal input voltage at the analog modulation input. The average power is defined as half the sum of maximum and minimum power.

Conditions: Modulation input signal as specified, modulation frequency as specified.

Measurement: Using a photoreceiver (of sufficient bandwidth) and an oscilloscope.

Frequency Response

Specifies the frequency $f_{\rm resp}$ up to which the power meter's transfer function is above a defined cutoff level.

The transfer function is the measured optical amplitude on a modulated sinusoidal optical input divided by the actual amplitude over modulation frequency, expressed in dB.




Figure 55 Frequency Response

Conditions: Cutoff level as specified. Other conditions as specified.

Insertion loss

Specifies the "Total loss" at a displayed attenuation of 0 dB

Conditions: Temperature range, wavelength range, source polarization input polarization, input mode (multimode attenuators only) as specified.

Measurement: With Fabry-Pérot Laser Source, polarization controller (if applicable) and power meter, using high quality connectors in perfect optical condition.

NOTE

"Insertion loss excluding connectors" means connector loss is not part of the loss value.

"Insertion loss including connectors" includes the additional insertion loss from the connectors.

Linearity

Specifies the relative difference between the measured (displayed) power ratio D_x/D_0 and the actual (true) power ratio P_x/P_0 , caused by changing the power level from the reference level P_0 to an arbitrary level P_x .

Powers expressed in Watts.

Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 74 and "Drift (dark)" (see, "Drift (dark)" on page 71).

$$L_{\%} = \left(\frac{D_x/D_0}{P_x/P_0} - 1\right) \text{ when expressed in \%.}$$

$$L_{\rm dB} = 10 \log \left(\frac{D_x / D_0}{P_x / P_0} \right)$$
 when expressed in dB.

Conditions: Reference power level $P_0 = 0.01 \text{ mW}$ (-20 dBm). Other conditions as specified.

Linewidth

The 3 dB width of the optical spectrum, expressed in Hertz.

Conditions: SBS suppression off.

Measurement:

Using a self-heterodyning technique: The output of the laser under test is sent through a Mach-Zehnder interferometer in which the length

Agilent N77xx Series

difference of the two arms is longer than the coherence length of the laser. The electrical noise spectrum of the photodetector current is measured with Agilent lightwave signal analyzer, and the linewidth is calculated from the heterodyne spectrum.

Alternatively, *Using a heterodyning technique:* The output of the laser under test is mixed with another laser of the same type on a wide bandwidth photodetector. The electrical noise spectrum of the photodetector current is measured with Agilent Lightwave signal analyzer, and the linewidth is calculated from the heterodyne spectrum. (Lightwave signal analyzer settings: resolution bandwidth 1 MHz, video bandwidth 10 kHz, sweep time 20 ms, single scan).

Maximum safe (input) power

The maximum input power that can be applied to any port of the power meter or attenuator without permanent change to its characteristics.

Attention: Applying more than the specified power might damage the power meter or attenuator!

	For input powers > 20 dBm, clean connectors in good condition are vital to avoid
ΝΟΤΕ	For input powers > 20 ubin, clean connectors in good condition are vital to avoid
NUL	thermally induced fiber damage.

Minimum averaging time

The minimum selectable averaging time.

Maximum output power

The maximum achievable output power of the tunable laser source and the maximum output power for which the tunable laser source specifications apply.

Conditions: As specified.

Measurement: Using a power meter at the output of the instrument.

Noise pp (dark)

Specifies the peak-to-peak (pp) change of displayed power over time with zero input power (dark), excluding "Drift (dark)" (see, "Drift (dark)" on page 71). Powers expressed in Watts.

Conditions: All ports dark. Averaging time as specified. Other conditions as specified.

Measurement: As for "Drift (dark)". From the measured Pi and the fitted curve F_i the noise curve $N_i = P_i - F_i$ is calculated.

The noise N is

NOTE

 $N = max_i\{N_i\} - min_i\{N_i\}$

Note: The noise pp specification (rather than noise 2 sigma) is given for legacy reasons.

Noise 2 sigma (dark)

Specifies twice the standard deviation of the measured dark power (with zero input power), excluding "Drift (dark)". Power are expressed in Watts.

Conditions: All ports dark. Averaging time as specified. Other conditions as specified.

Measurement: As for "Drift (dark)" (see, "Drift (dark)" on page 71). From the measurement result Pi and the fitted curve Fi the noise N is calculated by

 $N = 2^*StDev_i\{P_i - F_i\}$, where $StDev_i\{\}$ denotes the standard deviation over indices i.

Operating humidity

Humidity range where the instrument is designed to be operated.

The instrument must not be operated outside this range. If previously stored beyond this range, wait for acclimation before turning on the instrument.

Operating temperature

Temperature range for which the specifications apply if not differently stated.

The instrument must not be operated outside this range. If previously stored beyond this range, wait for acclimation before turning on the instrument.

Agilent N77xx Series

Operating temperature and humidity

The ambient temperature range and humidity range of the tunable laser source for which the specifications apply.

ΝΟΤΕ

If the tunable laser source module is rack- mounted the temperature and humidity within the rack apply.

Output isolation

The insertion loss of the built-in isolator in the backward direction.

Measurement: This characteristic cannot be measured from outside the module. It is based on known isolator characteristics.

Optical path blocking

Maximum possible attenuation.

Polarization dependent loss (PDL)

The dependence of the "Total loss" (see, "Total loss" on page 86) on the input polarization state, expressed as the full difference (in dB) between the highest and the lowest total loss.

Conditions: Temperature range, input mode as specified.

Measurement: Using high quality connectors in perfect optical condition.

NOTE

PDL is a peak-to-peak value.

Polarization dependent responsivity (PDR)

Specifies the dependence of the measured power on the input state of polarization (SOP) as \pm half the difference between the highest and the lowest measured power over all SOPs, expressed in dB.

Powers expressed in dBm.


Conditions: As specified.

Polarization extinction ratio

Specifies the ratio of the optical power in the slow axis of a connected polarization-maintaining fiber to optical power in the fast axis, expressed in dB

Conditions: Applicable to tunable laser sources utilizing polarization maintaining fiber that has its TE mode in the slow axis and aligned with the connector key.

Measurement: Using a polarization analyzer at the end of a polarization-maintaining patchcord, by sweeping the wavelength to create circular traces on the Poincaré sphere. Calculate the polarization extinction ratio from the diameters of these circles.



Port separation (static)

Specifies the difference between the power P applied to one port and the measured power D on any other port with no applied power (dark). Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 74) and "Drift (dark)" (see, "Drift (dark)" on page 71).

Powers expressed in dBm, port separation S expressed in dB,

S = P - D

Agilent N77xx Series

4 Specifications and Regulations Compliance




Figure 56 Port separation (static)

Conditions: Applied power level as specified.

Port separation, dynamic


Specifies the difference between step power input P_{step} (with positive or negative slope) applied to one port and the maximum measured power change D_{chng} (positive or negative) on any other port without applied power (dark). Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 74) and "Drift (dark)" (see, "Drift (dark)" on page 71). Powers expressed in dBm, dynamic port separation SD expressed in dB.

 $SD = P_{step} - D_{chng}$

Conditions: Power step height as specified (1 mW).

Power flatness versus wavelength

Specifies \pm half the span (in dB) between the maximum and the minimum actual power levels of the tunable laser source when changing the wavelength.



Conditions: Uninterrupted tunable laser source output power, constant power setting, constant temperature.

Power repeatability

The uncertainty in reproducing the power level after changing and resetting the power level. The power repeatability is ± half the span between the highest and lowest actual power (in dBm).


Conditions: Uninterrupted tunable laser source output power, constant wavelength setting, constant temperature.

NOTE

The long-term power repeatability can be obtained by taking the power repeatability and power stability into account.

Power stability

Specifies the change of the power level of the tunable laser source over time, expressed as ± half the span (in dB) between the highest and lowest actual power.



Conditions: Time span as specified. Uninterrupted tunable laser source output power, constant wavelength and power level settings, constant temperature.

Power range

Specifies the range from the smallest input power that causes a significant change of the measured (displayed) power to the highest power for which the specifications apply if not differently stated.

Measurement: The lower limit corresponds to half the "Noise pp (dark)" (see, "Noise pp (dark)" on page 74).

Agilent N77xx Series

Power setting mode

Operating mode where the user sets the desired output power. The power control function of the attenuator automatically adjusts the attenuation to get the set output power (see also, "Attenuation setting mode" on page 70).

Range in attenuation setting mode

Attenuation range that can be set at the attenuator in the Attenuation Setting Mode and for which the specifications apply (if not differently stated).

Range in power setting mode

Range of optical output power that can be set at the attenuator in the Power Setting Mode and for which the specifications apply (if not differently stated).

ΝΟΤΕ	The actual minimum value is determined by the optical input power minus the
NUL	maximum attenuation.

Range settling time (one range step)

When changing from one power meter range to the next (one range higher or lower), the settling time specifies the time required for the power meter to settle into the new range. Measurements starting after that period are within specifications.

The time starts when the instrument receives the range change command (manual range) or when the instrument detects the requirement for power meter range change (autorange).

Conditions: As specified.

Relative intensity noise (RIN)

Specifies the ratio between the mean-square of the optical power fluctuation amplitude Pf,B within a specified frequency range f and for bandwidth B, and the square of the average optical power P_{avg} .

$$RIN = \frac{\left\langle \Delta P_{f,B}^2 \right\rangle}{P_{avg}^2 \cdot B} \left[\frac{1}{Hz} \right]$$

RIN, if expressed as "dB/Hz", is calculated by:


$$RIN_{\rm dB/Hz} = 10 \log \left(\frac{\Delta P_{f,B}^2 \cdot 1 \,\rm{Hz}}{P_{\rm avg}^2 \cdot B} \right)$$


Conditions: As specified.

Measurement: Using an Agilent Lightwave signal analyzer and bandwidth set to 3 MHz.

Relative port to port uncertainty

Specifies \pm half the maximum difference between the measured powers of different ports when applying the same stimulus. Measured powers expressed in dBm, the relative port to port uncertainty *U* expressed in \pm dB.





Conditions: As specified.

Relative uncertainty of monitor power meter

When changing the output power of the attenuator, the relative power meter uncertainty is the maximum error of the displayed output power ratio to the actual output power ratio. This uncertainty is caused by the internal power meter's nonlinearity and noise and by errors in the sensing hardware, expressed as \pm half the span of all possible errors with an offset due to the noise level of the power meter. Symbol *RU*.

Conditions: Reference power level 1mW, wavelength and polarization state constant; power limitations as in the Specifications; zeroing prior to measurement.

Agilent N77xx Series

4 Specifications and Regulations Compliance

Measurement: The nonlinearity is calculated using the displayed power ratio D2 / D1 and the corresponding actual power ratio P2 / P1:

$$N_{dB}(1,2) = 10 \log \left(\frac{D_1 / D_2}{P_2 / P_1} \right)$$

The relative power meter uncertainty is then calculated using:


$$RU_{dB} = \pm \frac{max|N_{dB}(1,2)}{2}$$

RU is then expressed in dB and the offset is expressed in pW.

Absolute power accuracy attainable with the help of an external optical power meter.

Relative wavelength (frequency) accuracy

When randomly changing the wavelength (frequency) of the tunable laser source and measuring the differences between the displayed and the actual wavelength (frequency), the relative wavelength (frequency) accuracy is \pm half the span between the maximum and the minimum value of all differences.



Conditions: Uninterrupted tunable laser source output power, constant power setting, constant temperature.

Repeatability (of attenuation/total loss, or power level)

Specifies the uncertainty in reproducing the "Total loss" on page 86 (or "Attenuation" on page 69, or power level) after randomly changing

NOTE

www.valuetronics.com

82

and resetting the attenuation or power level. The repeatability is defined as \pm half the difference between the highest and the lowest measured attenuation or power level over the repeated changing and resetting of the attenuation or power level.




Figure 58 Repeatability

Conditions: Constant operating conditions. Operating in power setting mode requires zeroing prior to measurement.

Resolution

The smallest possible increment or decrement in Attenuation Setting Mode or Power Setting Mode.

Return loss

The ratio between the incident power and the power reflected back from the instrument into the input fiber, expressed in dB.

Conditions: As specified. (For attenuators: applicable to both ports with fiber cables with high quality connectors in perfect optical condition on both attenuator ports and the respective second port terminated (zero reflectance).)

Measurement: With return loss meter and non-coherent source of applicable wavelength.

NOTE

The measurement result includes internal reflections in the instrument, such as reflections from both attenuator ports.

SBS suppression - effective linewidth

Specifies the peak-to-peak change of the periodically modulated wavelength resulting from the SBS suppression feature, expressed in Hertz.

Side-mode suppression ratio

The ratio of optical power in the main mode to the optical power of the highest sidemode, expressed in dB:

$$SSR_{dB} = 10 \cdot \log\left(\frac{P_{signal}}{P_{highestsidemode}}\right)$$

Conditions: Within a distance from 0.1 to 6 GHz to the signal's optical frequency.


Measurement: Using the Agilent Lightwave signal analyzer, by analyzing the heterodyning between the main signal and the highest sidemode.


Settling Time

Maximum time needed to change the attenuation or power level by a specified step, from the beginning until the end of the change.

Signal to source spontaneous emission ratio

Specifies the ratio between signal power and maximum spontaneous emission (SSE) power. The SSE power is determined in a specified bandwidth within a ± 3 nm window around the signal wavelength, where ± 1 nm around the signal wavelength are excluded, expressed in dB per nm.





Measurement: Using an optical spectrum analyzer at 0.5 nm resolution bandwidth (to address the possibility of higher SSE within a narrower bandwidth), then extrapolated to 1 nm bandwidth.

Specification wavelength range

The wavelength range for which all specifications and characteristics apply, if not differently stated.

Spectral ripple (power meter; due to interference)

A coherent input signal causes optical interference1 between reflective interfaces within the power meter's optical assembly, including the end face of a connected fiber. Optical interference causes a wavelength-dependent periodic change of the power meter's responsivity. Spectral ripple is defined as \pm half the difference between the highest and the lowest power resulting from the periodic change.

Conditions: Constant input power, constant state of polarization (SOP). Connector as specified. Other conditions as specified.

Measurement: Using a laser source without spectral ripple or with known (monitored) spectral ripple and linewidth < 100 MHz and stepping the wavelength over the specified range. The periodic part in the measurement result is separated by subtracting the polynomial 2nd order fit (low frequency part).

Agilent N77xx Series

Storage conditions

Allowed temperature and humidity range for the non-operating instrument. Wait for acclimation to within the "Operating temperature" before turning on the instrument.

Total loss


The change of power level after inserting the attenuator between two connectorized patchcords, at an arbitrary attenuation setting, expressed in dB. Its value TL can be calculated from:

$$TL(dB) = 10\log \frac{P_a}{P_b} = P_a(dBm) - P_b(dBm)$$

where:

 P_a = power measured at the end of the two patchcords.

 $P_{\rm b}$ = power measured after the insertion of the attenuator.





Conditions: Patchcord cables with high quality connectors in perfect optical condition on both attenuator ports.

N O T E	The total loss depends on the attenuation setting.
NOTE	Total loss includes the loss of one additional connector pair.

Total uncertainty

The maximum relative difference U between measured (displayed) power M and actual (true) power P for a specified set of operating conditions, expressed in percent. Powers expressed in Watts. Excluding "Noise pp (dark)" (see, "Noise pp (dark)" on page 74) and "Drift (dark)" (see, "Drift (dark)" on page 71).

U = 100% * (M/P - 1)

Conditions: Operating conditions as specified. Other conditions as specified.

Transition speed

Specifies the settable attenuation change rate of the attenuator.

Tuning time

Specifies the time needed to tune to another wavelength.

Conditions: Tuning distance as specified. The time is measured from sending the command to the TLS until the tuning operation has finished.

Uncertainty at reference conditions

The maximum relative difference U between measured (displayed) power M and actual (true) power P for the specified set of reference conditions, expressed in percent. Powers expressed in Watts.

U = 100% * (M/P - 1)

Conditions: Reference conditions as specified.

Wavelength range

The range of wavelengths for which the power meter is calibrated, or that can be set at the attenuator and for which the specifications apply (if not differently stated).

Wavelength (frequency) repeatability

The random uncertainty in reproducing a wavelength (frequency) of the tunable laser source after changing and re-setting the wavelength (frequency). The wavelength (frequency) repeatability is ± half the span

Agilent N77xx Series

between the maximum and the minimum of all actual values of this wavelength (frequency).

Conditions: Uninterrupted tunable laser source output power, constant power setting, constant temperature.

ΝΟΤΕ

The long-term wavelength repeatability can be obtained by taking the wavelength repeatability and wavelength stability into account.

Wavelength resolution

The smallest selectable wavelength (frequency) increment or decrement.

Wavelength (frequency) stability

Specifies the change of the actual wavelength (frequency) of the tunable laser source over time, expressed as \pm half the span between the maximum and minimum of all wavelengths (frequencies).

Conditions: Time span as specified, uninterrupted tunable laser source output power, constant wavelength and power level settings, constant temperature.

Warm-up time

Time after power-up of the acclimated instrument after which the specifications and characteristics apply.

Literature

Fiber Optic Test and Measurement, Hewlett Packard Professional Books, edited by Prentice Hall, ISBN 0-13-534330-5

Technical Specifications - Optical Multiport Power Meter

	Agilent N7744A,	N7745A			
Sensor element	InGaAs	InGaAs			
Wavelength range	1250 nm to 1650	1250 nm to 1650 nm			
Specification wavelength range	1250 nm to 1625	nm (if not stat	ed differently)		
Power range	-80 dBm to +10 d	dBm			
Maximum safe power	+16 dBm				
Data logging capability	2 buffers per por points	2 buffers per port, each with capacitiy for 1 Mio. measurement points			
Averaging time	1 µs to 10 s				
Applicable fiber type	Standard SM and	Standard SM and MM \leq 62.5 μ m core size, NA \leq 0.24			
Uncertainty at reference conditions ^{[f1] [f3]}	±2.5%				
Total uncertainty ^{[f2] [f5] [f6]}	±4.5%				
Relative port to port uncertainty [f1] [f3] [f4] [f10]	typ. ±0.05 dB	typ. ±0.05 dB			
Linearity ^{[f5] [f6]} at (23 ± 5)°C	±0.02 dB ± 3 pW	±0.02 dB ± 3 pW			
over operating temperature	±0.04 dB ± 5 pW	±0.04 dB ± 5 pW			
	< ±0.015 dB (152	<pre></pre> <pre></pre> <pre></pre> <pre></pre> <pre></pre> <pre></pre> <pre></pre> <pre></pre> <pre></pre>			
Polarization dependent responsivity ^{[f3] [f7]}	$typ. < \pm 0.01 \text{ dB} (1250 \text{ nm to } 1580 \text{ nm})$				
		< ±0.01 dB (1520 nm to 1625 nm)			
Spectral ripple (due to interference) ^[f9]	typ. $< \pm 0.01$ dB (1250 nm to 1520 nm)				
Drift (dark) ^[f4]	±9 pW				
Noise peak-to-peak (dark) ^[f3]	·	aging time, 300	s observation tim	e)	
Noise 2σ ^[f3] (100,000 samples)	Averaging time:		25 µs	1 ms	
PM range -30 dBm	typ.	< 0.1 nW	< 0.025 nW	< 0.005 nW	
-20 dBm	typ.	< 1.5 nW	< 0.15 nW	< 0.02 nW	
-10 dBm	typ.	< 6 nW	< 0.5 nW	< 0.08 nW	
0 dBm	typ.	< 60 nW	< 4 nW	< 0.8 nW	
+10 dBm	typ.	< 600 nW	< 40 nW	< 8 nW	
Dynamic range (logging mode) ^{[f3] [f4]}	Averaging time:	1 µs	25 µs	1 ms	
PM range -30 dBm	typ.	> 43 dB	> 49 dB	> 57 dB	
-20 dBm	typ.	> 43 dB	> 54 dB	> 62 dB	
-10 dBm	typ.	> 46 dB	> 57 dB	> 64 dB	
0 dBm	typ.	> 46 dB	> 57 dB	> 63 dB	
+10 dBm	typ.	> 43 dB	> 54 dB	> 60 dB	
Port separation ^[f5]	> 85 dB (CW) (o	ne neighbor po	rt with 0 dBm)		
Port separation, dynamic ^[f5]	>70~dB typ. (one neighbor port with 0 dBm power in 0 dBm power meter range)				
Frequency response	3 dB cutoff frequ	iency at 1 µs av	veraging time, typ.		
-30 dBm range	10 kHz				
-20 dBm range	130 kHz				
-10 to +10 dBm range	250 kHz				

Agilent N77xx Series

4 Specifications and Regulations Compliance

	Agilent N7744A, N7745A		
Return loss [f8]	> 50 dB (1520nm to 1580nm) typ. > 57 dB (1280 nm to 1580 nm)		
Line Power	AC 100 - 240 V ± 10%, 50/60Hz, 60VA max.		
Operating temperature	+5 °C to +40 °C		
Operating humidity	15% to 95%, non-condensing		
Storage conditions	-40 °C to +70 °C		
Warm-up time	20 min.		
Recommended recalibration period	24 months		
Dimensions	372 mm × 212 mm × 43 mm (excluding front and back rubber cushions)		
Weight	3 kg (6 lb)		
LXI Compliance	LXI Class C, ver. 1.2		
 [f1] Reference conditions: Single mode fiber SMF 9 μm Power level: -20 dBm to 0 dBm On day of calibration (add ±0.3% for aging over year; add ±0.6% for aging over two years) Spectral width of source < 10 nm full width himum (FWHM) Wavelength setting of power sensor corresponsource wavelength ± 0.4 nm [f2] Operating conditions: Single mode fiber SMF. For multimode fiber, ty within one year of calibration; add ±0.3% for year Spectral width <10 nm FWHM 	alf maxi- [f6] Power range -60 dBm to +10 dBm [f7] Straight connector, SMF nds to [f8] Connector 8° angled, ceramic ferrule, SMF [f9] For constant state of polarization, source linewidtl <100 MHz, angled connector 8°, wavelength range 1520 nm to 1625 nm. Twical for 1250 nm to 1520 nm to 1625 nm.		
 Spectral width < 10 nm FWHW Wavelength setting of power sensor correspo source wavelength ± 0.4 nm 	nds to [f10] Same 4-detector block, same wavelength		

90

Declarations of Conformity - Optical Multiport Power Meter

		I	rding to EN ISO/IEC 17050-1:2004	
Manufacture Manufacture	r's Name: r's Address:	Agilent Techno No. 1 Yishun A SINGAPORE Singapore		nal) Pte. Ltd.
Declares un	der sole respoi	nsibility that the	product as originally deliv	ered
Product Nan Product Nun		N7744A	oort Power Meter 4 Channels 8 Channels	
Product Opt	ions:	This declaration	on covers all options of the a	bove system.
	th the essentia E marking acc		of the following applicable	European Directives, and
	•	tive 2006/95/EC		
	IC Directive 200			
and conform	is with the follo	wing product st	tandards:	
ЕМС	Standard	2005 / EN 61326-1:2	006	Limit
	IEC 61000 IEC 61000 IEC 61000 IEC 61000 IEC 61000 IEC 61000	:2003 / EN 55011:199)-4.2: 2001 / EN 6100)-4.3:2002+A1:2002 /)-4.4:2004 / EN 6100)-4.5: 2001 / EN 6100)-4.6:2003 / EN 6100)-4.6:2003 / EN 6100)0-4-5:1995+A1:2001 0-4-6:1996+A1:2001	Group 1 Class A 4 kV CD, 4 kV AD 3 V/m / 80 MHz-2 GHz, 1V/m / 2- 0.5 kV signal lines, 1 kV power lin 0.5 kV line-line, 1 kV line-ground 3 V, 0.15-80 MHz 0% for 1/0.5 (0°, 180°) cycle 0% for 250/300 cycles 70% for 25/30 cycles
		CES/NMB-001:2004 New Zealand: AS/NZS	S CISPR 11:2004	
Safety	IEC 61010-1:2001 / EN 61010-1:2001 Canada: CAN/CSA-C22.2 No. 61010-1:2004 USA: UL 61010-1:2004			
Supplement	ary Informatior	1:		
			ation with Agilent Technolog ced on the EU market after:	ies test systems.
2010-May-02			Hans-Martin Fisch	ner
Date			Name Product Quality &	Compliance

Year of first marking: 2008

Template:A5971-5302-2, Rev. F.00

Revision: B

Agilent N77xx Series

4 Specifications and Regulations Compliance

Regulations Information

EMC Canada

These ISM devices comply with Canadian ICES-001. Ces appareils ISM sont conformes à la norme NMB-001 du Canada.

Acoustic Noise Information Germany

ISO 7779 L_{pA} < 70 dB(A), normal operation, operator position.

Remove all doubt

Remove all doubt

Our repair and calibration services will get your equipment back to you, performing like new, when promised. You will get full value out of your Agilent equipment throughout its lifetime. Your equipment will be serviced by Agilent-trained technicians using the latest factory calibration procedures, automated repair diagnostics and genuine parts. You will always have the utmost confidence in your measurements. Agilent offers a wide range of additional expert test and measurement services for your equipment, including initial start-up assistance, onsite education and training, as well as design, system integration, and project management. For more information on repair and calibration services, go to

www.agilent.com/find/removealldoubt

Agilent E-mail Updates

Get the latest information on the products and applications you select.



www.agilent.com/find/emailupdates

Agilent Direct

Quickly choose and use your test equipment solutions with confidence.



www.agilent.com/find/agilentdirect

Agilent Open

Agilent Open simplifies the process of connecting and programming test systems to help engineers design, validate and manufacture electronic products. Agilent offers open connectivity for a broad range of systemready instruments, open industry software, PC-standard I/O and global support, which are combined to more easily integrate test system development.



www.agilent.com/find/open

LXI

LXI is the LAN-based successor to GPIB, providing faster, more efficient connectivity. Agilent is a founding member of the LXI consortium.



www.lxistandard.org

4 Specifications and Regulations Compliance

www.agilent.com

For more information on Agilent Technologies products, applications and services, please contact your local Agilent office. The complete list is available at:

www.agilent.com/find/contact us

Agilent N77xx Series



Accessories - Overview

The instruments are available in various configurations for the best possible match to the most common applications.

This chapter provides information on the available options and accessories.

Instrument and Options - Agilent N7744A and N7745A

N7744A, N7745A ordering options

N7744A	Optical Multiport Power Meter (4 channel)	
N7745A	Optical Multiport Power Meter (8 channel)	
includes: USB and cross-over l	LAN cables	
Accessories		
N7740FI	FC Connector Adapter for Optical Multiport Power Meter	
N7740KI	SC Connector Adapter for Optical Multiport Power Meter	
N7740LI	LC Connector Adapter for Optical Multiport Power Meter	
N7740MI	MU Connector Adapter for Optical Multiport Power Meter	
N7740ZI	Blank Zeroing Connector Adapter for Optical Multiport Power Meter	
N7740BI	Bare Fiber Connector Adapter for Optical Multiport Power Meter	
81004BM	Bare Fiber Holder Set for 0 - 400 μm Fibers (4 each with gauge)	
81009BM	Bare Fiber Holder Set for 400 - 900 μm Fibers (4 each with gauge)	
N7744A-100	Rack Mount Kit for one or two Multiport Power Meter Units	
Service and Repair R-51B-001-C	1 year Return-to-Agilent warranty	
R-51B-001-3C	1 year Return-to-Agilent warranty extended to 3 years	
R-51B-001-5C	1 year Return-to-Agilent warranty extended to 5 years	
R-50C-011-3	Agilent calibration up-front support plan 3 year coverage	
R-50C-011-5	Agilent calibration up-front support plan 5 year coverage	
R-50C-021-3	ANSI Z540-1-1994 up-front plan 3 year coverage	
	ANSI Z540-1-1994 up-front plan 5 year coverage	



Agilent Technologies

Cables

- GPIB Cable, 10833A, 1 m (3.3 feet)
- GPIB Cable, 10833B, 2 m (6.6 feet)
- GPIB Cable, 10833C, 4 m (13.2 feet)
- GPIB Cable, 10833D, 0.5 m (1.6 feet)
- GPIB Adapter, 10834A, 2.3 cm extender. Use this adapter if there is no space to connect your GPIB cable directly to a GPIB interface.
- USB Cable, 8121-1243, 2 m (6.6 feet)
- Cross-over LAN cable, 5061-0701, 2.1 m (7 feet)



Maintenance and Troubleshooting - Overview

This section provides detailed information on:

- Cleaning Instructions
- Firmware Upgrades
- Error Messages
- Data Sanitation

Cleaning Instructions

Cleaning Instructions

The following Cleaning Instructions contain some general safety precautions, which must be observed during all phases of cleaning. Consult your specific optical device manuals or guides for full information on safety matters.

Please try, whenever possible, to use physically contacting connectors, and dry connections. Clean the connectors, interfaces, and bushings carefully after use.

If you are unsure of the correct cleaning procedure for your optical device, we recommend that you first try cleaning a dummy or test device.

Agilent Technologies assume no liability for the customer's failure to comply with these requirements.



Cleaning Instructions for this Instrument

The Cleaning Instructions apply to a number of different types of Optical Equipment. Most relevant for this instrument is the section "How to clean connector adapters" on page 107 for cleaning the quadadapters, and "How to clean optical glass plates" on page 108 for cleaning the front of the instrument.

Safety Precautions

Please follow the following safety rules:

- Do not remove instrument covers when operating.
- Ensure that the instrument is switched off throughout the cleaning procedures.
- Use of controls or adjustments or performance of procedures other than those specified may result in hazardous radiation exposure.
- Make sure that you disable all sources when you are cleaning any optical interfaces.
- Under no circumstances look into the end of an optical device attached to optical outputs when the device is operational. The laser radiation is not visible to the human eye, but it can seriously damage your eyesight.
- To prevent electrical shock, disconnect the instrument from the mains before cleaning. Use a dry cloth, or one slightly dampened with water, to clean the external case parts. Do not attempt to clean internally.
- Do not install parts or perform any unauthorized modification to optical devices.
- Refer servicing only to qualified and authorized personnel.

Why is it important to clean optical devices?

In transmission links optical fiber cores are about 9 μ m (0.00035") in diameter. Dust and other particles, however, can range from tenths to hundredths of microns in diameter. Their comparative size means that they can cover a part of the end of a fiber core, and as a result will reduce the performance of your system.

Furthermore, the power density may burn dust into the fiber and cause additional damage (for example, 0 dBm optical power in a single mode fiber causes a power density of approximately 16 million W/m2). If this happens, measurements become inaccurate and non-repeatable.

Cleaning is, therefore, an essential yet difficult task. Unfortunately, when comparing most published cleaning recommendations, you will discover that they contain several inconsistencies. In this section, we want to suggest ways to help you clean your various optical devices, and thus significantly improve the accuracy and repeatability of your lightwave measurements.

What do I need for proper cleaning?

Some Standard Cleaning Equipment is necessary for cleaning your instrument. For certain cleaning procedures, you may also require certain Additional Cleaning Equipment.

Standard Cleaning Equipment

Before you can start your cleaning procedure you need the following standard equipment:

- Dust and shutter caps
- Isopropyl alcohol
- Cotton swabs
- Soft tissues
- Pipe cleaner
- Compressed air

Dust and shutter caps

All of Agilent Technologies' lightwave instruments are delivered with either laser shutter caps or dust caps on the lightwave adapter. Any cables come with covers to protect the cable ends from damage or contamination.

We suggest these protected coverings should be kept on the equipment at all times, except when your optical device is in use. Be careful when replacing dust caps after use. Do not press the bottom of the cap onto the fiber too hard, as any dust in the cap can scratch or pollute your fiber surface.

Agilent N77xx Series

If you need further dust caps, please contact your nearest Agilent Technologies sales office.

Isopropyl alcohol

This solvent is usually available from any local pharmaceutical supplier or chemist's shop.

If you use isopropyl alcohol to clean your optical device, do not immediately dry the surface with compressed air (except when you are cleaning very sensitive optical devices). This is because the dust and the dirt is solved and will leave behind filmy deposits after the alcohol is evaporated. You should therefore first remove the alcohol and the dust with a soft tissue, and then use compressed air to blow away any remaining filaments.

If possible avoid using denatured alcohol containing additives. Instead, apply alcohol used for medical purposes. Never try to drink this alcohol, as it may seriously damage to your health.

Do not use any other solvents, as some may damage plastic materials and claddings. Acetone, for example, will dissolve the epoxy used with fiber optic connectors. To avoid damage, only use isopropyl alcohol.

Cotton swabs

We recommend that you use swabs such as Q-tips or other cotton swabs normally available from local distributors of medical and hygiene products (for example, a supermarket or a chemist's shop). You may be able to obtain various sizes of swab. If this is the case, select the smallest size for your smallest devices.

Ensure that you use natural cotton swabs. Foam swabs will often leave behind filmy deposits after cleaning.

Use care when cleaning, and avoid pressing too hard onto your optical device with the swab. Too much pressure may scratch the surface, and could cause your device to become misaligned. It is advisable to rub gently over the surface using only a small circular movement.

Swabs should be used straight out of the packet, and never used twice. This is because dust and dirt in the atmosphere, or from a first cleaning, may collect on your swab and scratch the surface of your optical device.

Soft tissues

These are available from most stores and distributors of medical and hygiene products such as supermarkets or chemists' shops.

We recommend that you do not use normal cotton tissues, but multilayered soft tissues made from non-recycled cellulose. Cellulose tissues are very absorbent and softer. Consequently, they will not scratch the surface of your device over time.

Use care when cleaning, and avoid pressing on your optical device with the tissue. Pressing too hard may lead to scratches on the surface or misalignment of your device. Just rub gently over the surface using a small circular movement.

Use only clean, fresh soft tissues and never apply them twice. Any dust and dirt from the air which collects on your tissue, or which has gathered after initial cleaning, may scratch and pollute your optical device.

Pipe cleaner

Pipe cleaners can be purchased from tobacconists, and come in various shapes and sizes. The most suitable one to select for cleaning purposes has soft bristles, which will not produces scratches.

There are many different kinds of pipe cleaner available from tobacco shops.

The best way to use a pipe cleaner is to push it in and out of the device opening (for example, when cleaning an interface). While you are cleaning, you should slowly rotate the pipe cleaner.

Only use pipe cleaners on connector interfaces or on feed through adapters. Do not use them on optical head adapters, as the center of a pipe cleaner is hard metal and can damage the bottom of the adapter.

Your pipe cleaner should be new when you use it. If it has collected any dust or dirt, this can scratch or contaminate your device.

The tip and center of the pipe cleaner are made of metal. Avoid accidentally pressing these metal parts against the inside of the device, as this can cause scratches.

Compressed air

Compressed air can be purchased from any laboratory supplier.

It is essential that your compressed air is free of dust, water and oil. Only use clean, dry air. If not, this can lead to filmy deposits or scratches on the surface of your connector. This will reduce the performance of your transmission system.

When spraying compressed air, hold the can upright. If the can is held at a slant, propellant could escape and dirty your optical device. First spray into the air, as the initial stream of compressed air could contain

Agilent N77xx Series

6 Maintenance and Troubleshooting

some condensation or propellant. Such condensation leaves behind a filmy deposit.

Please be friendly to your environment and use a CFC-free aerosol.

Additional Cleaning Equipment

Some Cleaning Procedures need the following equipment, which is not required to clean each instrument:

- Microscope with a magnification range about 50X up to 300X
- Ultrasonic bath
- Warm water and liquid soap
- Premoistened cleaning wipes
- Polymer film
- Infrared Sensor Card

Microscope with a magnification range about 50X up to 300X

A microscope can be found in most photography stores, or can be obtained through or specialist mail order companies. Special fiberscopes are available from suppliers of splicing equipment.

Ideally, the light source on your microscope should be very flexible. This will allow you to examine your device closely and from different angles.

A microscope helps you to estimate the type and degree of dirt on your device. You can use a microscope to choose an appropriate cleaning method, and then to examine the results. You can also use your microscope to judge whether your optical device (such as a connector) is severely scratched and is, therefore, causing inaccurate measurements.

Ultrasonic bath

Ultrasonic baths are also available from photography or laboratory suppliers or specialist mail order companies.

An ultrasonic bath will gently remove fat and other stubborn dirt from your optical devices. This helps increase the life span of the optical devices.

Only use isopropyl alcohol in your ultrasonic bath, as other solvents may damage.

Warm water and liquid soap

Only use water if you are sure that there is no other way of cleaning your optical device without corrosion or damage. Do not use hot water, as this may cause mechanical stress, which can damage your optical device.

Ensure that your liquid soap has no abrasive properties or perfume in it. You should also avoid normal washing-up liquid, as it can cover your device in an iridescent film after it has been air-dried.

Some lenses and mirrors also have a special coating, which may be sensitive to mechanical stress, or to fat and liquids. For this reason we recommend you do not touch them.

If you are not sure how sensitive your device is to cleaning, please contact the manufacturer or your sales distributor.

Premoistened cleaning wipes

Use pre-moistened cleaning wipes as described in each individual cleaning procedure. Cleaning wipes may be used in every instance where a moistened soft tissue or cotton swab is applied.

Polymer film

Polymer film is available from laboratory suppliers or specialist mail order companies.

Using polymer film is a gentle method of cleaning extremely sensitive devices, such as reference reflectors and mirrors.

Infrared Sensor Card

Infrared sensor cards are available from laboratory suppliers or specialist mail order companies.

With this card you are able to control the shape of laser light emitted. The invisible laser beam is projected onto the sensor card, then becomes visible to the normal eye as a round spot.

Take care never to look into the end of a fiber or any other optical component, when they are in use. This is because the laser can seriously damage your eyes.

Agilent N77xx Series

Preserving Connectors

Listed below are some hints on how best to keep your connectors in the best possible condition.

Making Connections

Before you make any connection you must ensure that all cables and connectors are clean. If they are dirty, use the appropriate cleaning procedure.

When inserting the ferrule of a patchcord into a connector or an adapter, make sure that the fiber end does not touch the outside of the mating connector or adapter. Otherwise you will rub the fiber end against an unsuitable surface, producing scratches and dirt deposits on the surface of your fiber.

Dust Caps and Shutter Caps

Be careful when replacing dust caps after use. Do not press the bottom of the cap onto the fiber as any dust in the cap can scratch or dirty your fiber surface.

When you have finished cleaning, put the dust cap back on, or close the shutter cap if the equipment is not going to be used immediately.

Keep the caps on the equipment always when it is not in use.

All of Agilent Technologies' lightwave instruments and accessories are shipped with either laser shutter caps or dust caps. If you need additional or replacement dust caps, contact your nearest Agilent Technologies Sales/Service Office.

Immersion Oil and Other Index Matching Compounds

Where it is possible, do not use immersion oil or other index matching compounds with your device. They are liable to impair and dirty the surface of the device. In addition, the characteristics of your device can be changed and your measurement results affected.

Cleaning Instrument Housings

Use a dry and very soft cotton tissue to clean the instrument housing and the keypad. Do not open the instruments as there is a danger of electric shock, or electrostatic discharge. Opening the instrument can

cause damage to sensitive components, and in addition your warranty will be voided.

Never open the instruments as they can be damaged. Opening the instruments puts you in danger of receiving an electrical shock from your device, and renders your warranty void.

Which Cleaning Procedure should I use?

Light dirt

If you just want to clean away light dirt, observe the following procedure for all devices:

- 1 Use compressed air to blow away large particles.
- **2** Clean the device with a dry cotton swab.
- **3** Use compressed air to blow away any remaining filament left by the swab.

Heavy dirt

If the above procedure is not enough to clean your instrument, follow one of the procedures below. Please consult "Cleaning Instructions for this Instrument" on page 98 for the procedure relevant for this instrument.

If you are unsure of how sensitive your device is to cleaning, please contact the manufacturer or your sales distributor

How to clean connectors

Cleaning connectors is difficult as the core diameter of a single-mode fiber is only about 9 μ m. This generally means you cannot see streaks or scratches on the surface. To be certain of the condition of the surface of your connector and to check it after cleaning, you need a microscope.

In the case of scratches, or of dust that has been burnt onto the surface of the connector, you may have no option but to polish the connector. This depends on the degree of dirtiness, or the depth of the scratches. This is a difficult procedure and should only be performed by skilled personal, and as a last resort as it wears out your connector.

Agilent N77xx Series

WARNING

Never look into the end of an optical cable that is connected to an active source.

To assess the projection of the emitted light beam you can use an infrared sensor card. Hold the card approximately 5 cm from the output of the connector. The invisible emitted light is project onto the card and becomes visible as a small circular spot.

Preferred Procedure

Use the following procedure on most occasions.

- 1 Clean the connector by rubbing a new, dry cotton-swab over the surface using a small circular movement.
- **2** Blow away any remaining lint with compressed air.

Procedure for Stubborn Dirt

Use this procedure particularly when there is greasy dirt on the connector:

- 1 Moisten a new cotton-swab with isopropyl alcohol.
- 2 Clean the connector by rubbing the cotton-swab over the surface using a small circular movement.
- **3** Take a new, dry soft-tissue and remove the alcohol, dissolved sediment and dust, by rubbing gently over the surface using a small circular movement.
- **4** Blow away any remaining lint with compressed air.

An Alternative Procedure

A better, more gentle, but more expensive cleaning procedure is to use an ultrasonic bath with isopropyl alcohol.

- **1** Hold the tip of the connector in the bath for at least three minutes.
- **2** Take a new, dry soft-tissue and remove the alcohol, dissolved sediment and dust, by rubbing gently over the surface using a small circular movement.
- **3** Blow away any remaining lint with compressed air.

How to clean connector adapters

Preferred Procedure

Use the following procedure on most occasions.

- 1 Clean the adapter by rubbing a new, dry cotton-swab over the surface using a small circular movement.
- 2 Blow away any remaining lint with compressed air.

Procedure for Stubborn Dirt

Use this procedure particularly when there is greasy dirt on the adapter:

- **1** Moisten a new cotton-swab with isopropyl alcohol.
- 2 Clean the adapter by rubbing the cotton-swab over the surface using a small circular movement.
- **3** Take a new, dry soft-tissue and remove the alcohol, dissolved sediment and dust, by rubbing gently over the surface using a small circular movement.
- 4 Blow away any remaining lint with compressed air.

How to clean bare fiber adapters

Bare fiber adapters are difficult to clean.

Protect from dust unless they are in use.

CAUTION

Never use any kind of solvent when cleaning a bare fiber adapter as solvents can damage the foam inside some adapters.

They can deposit dissolved dirt in the groove, which can then dirty the surface of an inserted fiber.

Preferred Procedure

Use the following procedure on most occasions.

1 Blow away any dust or dirt with compressed air.

How to clean optical glass plates

The multiport power meter has an optical glass plate to protect the sensor. This has a special coating that is sensitive to solvents, grease, liquid and mechanical abrasion. Take extra care when cleaning this plate.

Preferred Procedure

Use the following procedure on most occasions.

- 1 Clean the plate by rubbing a new, dry cotton-swab over the surface using a small circular movement.
- 2 Blow away any remaining lint with compressed air.

Procedure for Stubborn Dirt

Use this procedure particularly when there is greasy dirt on the plate:

- **1** Moisten a new cotton-swab with isopropyl alcohol.
- **2** Clean the plate by rubbing the cotton-swab over the surface using a small circular movement.
- **3** Using a new, dry cotton-swab remove the alcohol, any dissolved sediment and dust.
- 4 Blow away any remaining lint with compressed air.

Additional Cleaning Information

The following cleaning procedures may be used with other optical equipment:

How to clean bare fiber ends

How to clean bare fiber ends

Bare fiber ends are often used for splices or, together with other optical components, to create a parallel beam. The end of a fiber can often be scratched. You make a new cleave. To do this:

- **1** Strip off the cladding.
- 2 Take a new soft-tissue and moisten it with isopropyl alcohol.
- **3** Carefully clean the bare fiber with this tissue.

4 Make your cleave and immediately insert the fiber into your bare fiber adapter in order to protect the surface from dirt.

Other Cleaning Hints

Selecting the correct cleaning method is an important element in maintaining your equipment and saving you time and money. This Appendix highlights the main cleaning methods, but cannot address every individual circumstance.

This section contain some additional hints which we hope will help you further. For further information, please contact your local Agilent Technologies representative.

Making the connection

Before you make any connection you must ensure that all lightwave cables and connectors are clean. If not, then use appropriate the cleaning methods.

When you insert the ferrule of a patchcord into a connector or an adapter, ensure that the fiber end does not touch the outside of the mating connector or adapter. Otherwise, the fiber end will rub up against something which could scratch it and leave deposits.

Lens cleaning papers

Note that some special lens cleaning papers are not suitable for cleaning optical devices like connectors, interfaces, lenses, mirrors and so on. To be absolutely certain that a cleaning paper is applicable, please ask the salesperson or the manufacturer.

Firmware Upgrades

Updated versions of the firmware for your instrument may be released from time to time. These will then be available from the Agilent web page (www.agilent.com/find/octfirmware) and will include the update application and instructions.

Agilent N77xx Series

Error Messages

SYST:ERR?

The syst:err? command returns the next error from the error queue (refer to the Programming Guide for details).

Each error has an error code and a short description of the error, separated by comma. Please note that the error codes are returned as signed integer numbers in a range from -32768 to +32767 (INT16). Negative error numbers are defined by the SCPI standard, positive error numbers are device dependent.

Data Sanitation

Memory Inventory

The instruments have the following memory size:

EEPROM 1Mb + 1Mb for each block of ports That is a N7745A, 2-port power meter has 3Mb. This memory does not contain any customer data.

Flash ROM 16MB This memory can contain calibration data, and so on, as well as user parameter settings, zeroing data and logged information (the log contains limited operational information).

DRAM volatile working memory. This contains data for current and recent calibrating and measuring.

Clearing User Data

Clearing the DRAM

All data in the DRAM is cleared by powering the instrument off.

Clearing the Flash ROM

The user data can be cleared from the Flash ROM by:

- **1** Reset the instrument to factory default.
 - **a** Send the command:

SYSTem:PRESet

OR

- **b** Press the LAN Reset button on the front panel briefly (<1 s).
- **2** Clear the instrument log.
 - **a** Send the special command

DIAG:SYST:FLAS:LOGF:ERAS

The contents of the log can be checked with the query READ:LOGF?.

- **3** Clear the Ethernet settings
 - **a** Connect to the instrument by some means other than the LAN, using VISA Assistant (for example over USB). Send the command:

SYSTem: COMMunicate: ETHernet: RESet

OR

b With the instrument connected over the LAN (the LAN LED is green), press and hold the LAN Reset button for more than 3s.

In either case, the instrument reboots.

Agilent N77xx Series

6 Maintenance and Troubleshooting

Agilent N77xx Series

Index

A

Absolute wavelength (frequency) accuracy, 69 Accessories, 95 Adapters, 33 Additional Cleaning Equipment, 99 Agilent N7744A, 7 Agilent N7745A, 7 Agilent N7751A, 8 Agilent N7752A, 8 Agilent N7766A, 9 Agilent N7768A, 9 Angled contact connector, 25 Attenuation, 69 Attenuation accuracy (uncertainty), 70 Attenuation setting mode, 70 Attenuator. 55 Auto Gain, 60 Averaging time, 70 averaging time, 60

B

Bare fiber adapters, 107 BNC connectors, 28, 62

C

Cables, 33 Cleaning connectors, 105 Cleaning Instructions, 97 Cleaning Procedures, 102 Compact multichannel attenuators, 8 Compact Tunable Laser, 10 Connect angled contact fiber end connectors, 25 Connector, 33 connectors, 104 Constant operating conditions, 71 Coupler, 55

D

Data logging capability, 71 DRAM, 110 Drift (dark), 71 Dynamic range (logging mode), 71

E

EEPROM, 110 Electrical Connectors, 7, 28 Electrical offset, 57 External analog modulation modulation depth, 72

F

Flash ROM, 110 Frequency Response, 72 Front Panel Indicators, 35 Front Panels - Multiport Power Meters, 21 Front Panels - Optical Attenuators, 21, 22 Front Panels - Tunable Laser Source, 21, 24

G

GPIB Interface, 7

GPIB interface, 33

Infrared sensor cards, 103 Initial Inspection, 7 Input trigger functions, 63 Input/Output Signals, 15 Insertion loss, 73 Instrument Markings, 16 IP address, 29 Isopropyl alcohol, 100

L

LAN Interface, 7 Laser class 1M label, 14 laser safety labels, 14 Line Power Cable, 17 Line Power Connectors, 15 Line Power Requirements, 13 Linearity, 73 Linewidth, 73 Logging, 65

Μ

Maximum output power, 74 Maximum safe (input) power, 74 Microscope, 102 Minimum averaging time, 74 MinMax, 60 Multimode attenuators, 25 Multiport Power Meter, 63

Agilent N77xx Series

Ν

N77xx Viewer, 56 Noise 2 sigma (dark), 75 Noise pp (dark), 74

0

Operating Environment, 13 Operating humidity, 75 Operating temperature, 75 Operating temperature and humidity, 76 Optical connectors. 25 Optical attenuators, 25 98 Optical Equipment, optical glass plate, 108 Optical Multiport Power Meter, 62 Optical path blocking, 76 Optical Power Meters, 57 Options, 95 Output isolation, 76 Output trigger functions, 63

P

Pipe cleaners, 101 Polarization Dependent Loss, 61 Polarization dependent loss (PDL), 76 Polarization extinction ratio, 77 Polymer film, 103 Port separation (static), 77 Port separation, dynamic, 78 Power flatness versus wavelength, 78 Power Meter, 57 Power Meters, 8 Power range, 79 Power repeatability, 79 Power setting mode, 80 Power stability, 79

Power Unit, 54 Pre-moistened cleaning, 103

R

Range in attenuation setting mode, 80 Range in power setting mode, 80 Range setting, 59 Range settling time, 80 ranging modes, 58 Rear Panel, 21, 24 Relative intensity noise, 80 Relative port to port uncertainty, 81 Relative uncertainty of monitor power meter, 81 Relative wavelength (frequency) accuracy, 82 Repeatability (of attenuation/total loss, or power level), 82 Resolution, 83 Return loss. 83

S

Safety, 11 Safety Considerations, 7, 11 Safety Information for Tunable Laser Modules, 13 Safety Symbols, 11 SBS suppression - effective linewidth, 84 Settling Time, 84 Side-mode suppression ratio, 84 Signal to source spontaneous emission ratio, 84 Specification wavelength range, 85 Spectral ripple (power meter; due to interference), 85 Stability, 65 Standard Cleaning Equipment, 99

Storage conditions, 86 Straight contact connectors, 25 Switch, 55

T

Technical Specifications - Optical Multiport Power Meter, 89 Temperature for operation and storage, 19 Total loss, 86 Total uncertainty, 87 Transition speed, 87 Trigger, 62 Tunable Laser Source, 69 Tuning time, 87

U

Ultrasonic baths, 102 Uncertainty at reference conditions, 87 USB Interface, 7

W

Warm-up time, 88 Watts, 54 Wavelength (frequency) repeatability, 87 Wavelength (frequency) stability, 88 Wavelength range, 87 Wavelength resolution, 88 Web Control Interface, 61

Ζ

Zeroing, 57

www.agilent.com

 $\ensuremath{\mathbb{C}}$ Agilent Technologies Deutschland GmbH 2010 Fourth edition, July 2010



N7744-90B01

